

Narrow Mouth Bottle, PP PP 세구병

Cat. No.	Model	용량	Size ($\varnothing \times h$, mm)	Neck id (mm)	PK Price ₩		CS Price ₩	
B08-46-007	T582070	4 ml	16×41	8.3	12ea	3,000	6pk	17,200
B08-46-014	T582080	8 ml	25×44	13	12ea	3,400	6pk	19,300
B08-46-021	T582090	15 ml	25×58	13.3	12ea	5,300	6pk	30,600
B08-46-028	T582100	30 ml	36×62	13.5	12ea	6,700	6pk	38,400
B08-46-035	T582110	60 ml	40.5×85	14.5	12ea	8,900	6pk	50,800
B08-46-042	T582120	125 ml	51×99	17.5	12ea	10,300	6pk	59,000
B08-46-049	T582130	250 ml	61×133		12ea	15,900	6pk	90,600
B08-46-056	T582140	500 ml	73×170	21.3	12ea	22,500	4pk	85,700
B08-46-063	T582150	1 L	92×216	28	6ea	21,400	4pk	81,400
B08-46-066	T582180	2 L	118×252	27.5	1ea	11,800	6ea	67,200

- 재질 : PP, -10 ~ 135°C 내열, Autoclavable
- USP Class VI에 따른 Medical Grade Resin 사용
- Integral Seal Ring이 있는 Linerless Leak Proof Closure
- 반투명 재질로 화학적 내성이 우수함
- Tarsons



HDPE Narrow Mouth Bottle HDPE 세구병

Cat. No.	Model	용량	Size ($\varnothing \times h$, mm)	Neck id (mm)	PK Price ₩		CS Price ₩	
B08-46-707	T583170	4 ml	16×41	8.3	12ea	3,000	6pk	17,200
B08-46-714	T583180	8 ml	25×44	13	12ea	3,400	6pk	19,300
B08-46-721	T583190	15 ml	25×58	13.3	12ea	5,800	6pk	33,000
B08-46-728	T583100	30 ml	36×62	13.5	12ea	7,000	6pk	39,900
B08-46-735	T583110	60 ml	40.5×85	14.5	12ea	8,800	6pk	50,100
B08-46-742	T583120	125 ml	51×99	17.5	12ea	10,300	6pk	59,000
B08-46-749	T583130	250 ml	61×133		12ea	15,900	6pk	90,600
B08-46-756	T583140	500 ml	73×170	21.3	12ea	22,500	4pk	85,700
B08-46-763	T583150	1 L	92×216	28	6ea	21,400	4pk	81,400

- 재질 : HDPE, -100 ~ 120°C 내열, Non-autoclavable
- USP Class VI에 따른 Medical Grade Resin 사용
- Integral Seal Ring이 있는 Linerless Leak Proof Closure
- 불투명 HDPE로 저온 보관용으로 알맞음
- Tarsons



Amber Narrow-Mouth Bottle, HDPE 갈색 HDPE 세구병

Cat. No.	Model	용량	Size ($\varnothing \times h$, mm)	Neck id (mm)	PK Price ₩		CS Price ₩	
B08-46-070	T581170	4 ml	16×41	8.3	12ea	3,600	6pk	20,600
B08-46-077	T581180	8 ml	25×44	13	12ea	4,400	6pk	25,400
B08-46-084	T581190	15 ml	25×58	13.3	12ea	5,700	6pk	32,300
B08-46-091	T581200	30 ml	36×62	13.5	12ea	8,400	6pk	48,100
B08-46-098	T581210	60 ml	40.5×85	14.5	12ea	10,100	6pk	57,700
B08-46-105	T581220	125 ml	51×99	17.5	12ea	13,200	6pk	75,500
B08-46-112	T581230	250 ml	61×133		12ea	20,000	6pk	114,000
B08-46-119	T581240	500 ml	73×170	21.3	12ea	28,800	4pk	109,900
B08-46-126	T581250	1 L	92×216	28	6ea	27,200	4pk	103,500

- 재질 : HDPE, -100 ~ 120°C 내열, Non-autoclavable
- USP Class VI에 따른 Medical Grade Resin 사용
- Integral Seal Ring이 있는 Linerless Leak Proof Closure
- 완전 불투명 HDPE재질로 빛에 민감한 제품의 보관용으로 알맞음
- Comply with USP 23.661 Requirement for Max. Light Transmission
- Tarsons

