

Wide Mouth Bottle, PP

PP 광구병

Cat. No.	Model	용량	Size (Ø×h, mm)	Neck id (mm)	PK Price ₩		CS Price ₩	
B08-46-133	T582200	30 ml	36.5×62.5	21.5	12ea	7,400	6pk	42,500
B08-46-140	T582210	60 ml	38.5×85		12ea	9,000	6pk	51,600
B08-46-147	T582220	125 ml	50.5×99	28.5	12ea	12,400	6pk	70,700
B08-46-154	T582230	250 ml	61.5×133	33	12ea	19,400	6pk	110,600
B08-46-161	T582240	500 ml	72.5×170	43.5	12ea	27,600	4pk	105,300
B08-46-168	T582250	1 L	91.5×199	53	6ea	27,600	4pk	105,300
B08-46-175	T582260	2 L	119×245	90	1ea	10,100	6pk	57,500

- 재질 : PP, -10 ~ 135°C 내열, Autoclavable
- USP Class VI에 따른 Medical Grade Resin 사용
- Integral Seal Ring이 있는 Linerless Leak Proof Closure
- 반투명 재질로 화학적 내성이 우수함 • Tarsons



HDPE Wide Mouth Bottle

HDPE 광구병

Cat. No.	Model	용량	Size (Ø×h, mm)	Neck id (mm)	PK Price ₩		CS Price ₩	
B08-46-833	T584200	30 ml	36.5×62.5	21.5	12ea	7,400	6pk	42,500
B08-46-840	T584210	60 ml	38.5×85		12ea	9,000	6pk	51,600
B08-46-847	T584220	125 ml	50.5×99	28.5	12ea	12,400	6pk	70,700
B08-46-854	T584230	250 ml	61.5×133	33	12ea	19,400	6pk	110,600
B08-46-861	T584240	500 ml	72.5×170	43.5	12ea	27,600	4pk	105,300
B08-46-868	T584250	1 L	91.5×199	53	6ea	27,900	4pk	106,300
B08-46-875	T584260	2 L	119×245	90	1ea	10,100	6pk	57,500

- 재질 : HDPE, -100 ~ 120°C 내열, Non-autoclavable
- USP Class VI에 따른 Medical Grade Resin 사용
- Integral Seal Ring이 있는 Linerless Leak Proof Closure
- 불투명 HDPE로 저온 보관용으로 알맞음 • Tarsons



Amber Wide Mouth Bottle, HDPE

갈색 HDPE 광구병

Cat. No.	Model	용량	Size (Ø×h, mm)	Neck id (mm)	PK Price ₩		CS Price ₩	
B08-46-182	T581300	30 ml	36.5×62.5	21.5	12ea	9,900	6pk	56,300
B08-46-189	T581310	60 ml	38.5×85		12ea	11,700	6pk	66,600
B08-46-196	T581320	125 ml	50.5×99	28.5	12ea	14,600	6pk	83,700
B08-46-203	T581330	250 ml	61.5×133	33	12ea	22,600	6pk	129,200
B08-46-210	T581340	500 ml	72.5×170	43.5	12ea	33,700	4pk	128,200
B08-46-217	T581350	1 L	91.5×199	53	6ea	38,500	4pk	146,500
B08-46-224	T581360	2 L	119×245	90	1ea	12,900	6pk	73,500

- 재질 : HDPE, -100 ~ 120°C 내열, Non-autoclavable
- USP Class VI에 따른 Medical Grade Resin 사용
- Integral Seal Ring이 있는 Linerless Leak Proof Closure
- 완전 불투명 HDPE 재질로 빛에 민감한 제품의 보관용으로 알맞음
- Comply with USP 23.661 Requirement for Max. Light Transmission
- Tarsons



Wide Mouth Square Bottle, PP

PP 사각 광구병

Cat. No.	Model	용량	Size (Ø×h, mm)	Neck id	PK Price ₩		CS Price ₩	
B08-46-470	T582221	60 ml	35.4×81.5	21 mm	12ea	15,500	6pk	88,400
B08-46-473	T582222	175 ml	52.1×100.7	28.7 mm	12ea	22,800	6pk	130,400
B08-46-476	T582223	250ml	60.3×112.4	33 mm	12ea	29,200	6pk	167,200
B08-46-479	T582224	500ml	75×141.1	43.8 mm	12ea	40,700	4pk	155,100
B08-46-482	T582225	1 L	90.7×178	53.2 mm	6ea	38,700	4pk	147,500

- 재질 : PP, -10 ~ 135°C 내열, Autoclavable
- USP Class VI에 따른 Medical Grade Resin 사용
- Integral Seal Ring이 있는 Linerless Leak Proof Closure
- 반투명 재질로 화학적 내성이 우수함 • Tarsons

