

Centrifuge Bottle, PP 원심 분리기용 PP 광구병

Cat. No.	Model	용량	Size (Ø×h)	RCF	PK Price ₩
B08-06-005	T544010	250 ml	59.6×127 mm	13200 g	4ea 21,600
B08-06-010	T544020	500 ml	69.3×166 mm	13700 g	4ea 26,300
B08-06-015	T544030	1 L	96.7×173 mm	7100 g	4ea 56,400
B08-06-020	T544040	1 L	96.7×167 mm		4ea 56,400

- 재질 : PP, -10 ~ 135°C 내열, Autoclavable
- 반투명 재질로 화학적 내성이 우수함
- 최적의 결과를 얻기 위하여 사용시 최저 80% 이상 Filling 되어야 함
- Tarsons



- Centrifuge Tube는 1127 ~ 1132 Page 참조

Vacuum Bottle, PP PP 진공병

Cat. No.	Model	용량	Size (Ø×h)	Cap	EA Price ₩
B08-46-245	T583255	2 L	118×270 mm	53 mm	16,600
B08-46-247	T582255	4 L	155×338 mm	83 mm	30,500

- 재질 : PP, -10 ~ 135°C 내열, Autoclavable
- USP Class VI에 따른 Medical Grade Resin 사용
- Heavy Duty Type으로 Thick-Wall을 가지고 있음, TPE Gasket
- Translucent, Leakproof • Neck id : 38.5 mm (2 L), 65 mm (4 L)
- Full Vacuum에서 24시간 동안 유지 할 수 있음 (20°C)
- Tarsons

Optional Accessory

Cat. No.	Model	Cap Size	EA Price ₩
B08-46-252	T583216	53 mm	39,200
B08-46-253	T583217	83 mm	62,500

- 재질 : PP Body, Silicon Tubing & Cap
- Port 수 : 3개



Optional Accessory

- Quick Filling & Venting Closure는 160 Page 참조



Clean Bottle, PE PE 크린 바틀

- 초순수로 세척하고, 크린 포장한 제품임
- Class 100 Cleanroom용
- Universal

세구 크린 바틀

Cat. No.	Model	용량	Size	PK Price ₩
B08-24-601	16260	100 ml	Ø48×h97 mm	5ea 31,500
B08-24-605	16261	250 ml	Ø63×h127 mm	5ea 32,400
B08-24-609	16262	500 ml	Ø78×h167 mm	2ea 19,800
B08-24-613	16263	1 L	Ø99×h202 mm	1ea 18,600

광구 크린 바틀

Cat. No.	Model	용량	Size	PK Price ₩
B08-24-666	6677	100 ml	Ø51×h90 mm	5ea 23,300
B08-24-670	6676	250 ml	Ø63×h121 mm	5ea 33,000
B08-24-674	6678	500 ml	Ø78×h151 mm	2ea 18,000
B08-24-678	6679	1 L	Ø100×h186 mm	1ea 15,800

