

Straight Side Round Bottle / Jar, PP
PP 단형 대 광구병 / 자

Cat. No.	Model	용량	Size	PK Price ₩		CS Price ₩	
B08-134-495	PLC03693	30 ml	Ø39×h42 mm	12ea	28,000	6pk	157,000
B08-134-500	PLC03697	60 ml	Ø49×h48 mm	12ea	29,500	4pk	110,100
B08-134-505	PLC03701	120 ml	Ø53×h72 mm	12ea	42,300	3pk	118,400
B08-134-510	PLC03704	250 ml	Ø84×h62 mm	12ea	69,600	3pk	195,100
B08-134-515	PLC03717	500 ml	Ø84×h99 mm	6ea	42,100	4pk	157,200
B08-134-520	PLC03720	1 L	Ø115×h99 mm	6ea	71,600	4pk	267,600

- 재질 : PP, -10 ~ 135°C 내열, Autoclavable
- Qorpak



Straight Side wide Mouth Jar / Bottle, PP
PP 단형 대 광구병 / 자

Cat. No.	Model	용량	Size (Ø×h, mm)	Neck id (mm)	PK Price ₩		CS Price ₩	
B08-46-601	T882080	125 ml	62×63	70	4ea	9,000	6pk	51,600
B08-46-603	T882090	250 ml	62×103		4ea	11,700	6pk	66,600
B08-46-605	T882200	500 ml	110×76	120	4ea	19,800	4pk	75,500
B08-46-607	T882222	1 L	110×135		4ea	27,200	4pk	103,500

- 재질 : PP, -10 ~ 135°C 내열, Autoclavable
- 반투명 재질로 화학적 내성이 뛰어남
- Tarsons



Straight Side wide Mouth Jar / Bottle, PMP
PMP 단형 대 광구병 / 자

Cat. No.	Model	용량	Size (Ø×h, mm)	Neck id (mm)	PK Price ₩		CS Price ₩	
B08-46-611	T883080	125 ml	62×63	70	4ea	18,600	6pk	104,300
B08-46-613	T883090	250 ml	62×103		4ea	24,200	6pk	135,700
B08-46-615	T883300	500 ml	110×76	120	4ea	47,100	4pk	175,900
B08-46-617	T883333	1 L	110×135		4ea	59,200	4pk	221,200

- 재질 : PMP Body, 0 ~ 150°C 내열, Autoclavable / PP Cover
- 투명 재질로 화학적 내성이 뛰어남
- Museum에서 사용에 알맞음(Formaldehyde loss is only 2% per year)
- Tarsons



Straight Side wide Mouth Jar / Bottle, PC
PC 단형 대 광구병 / 자

Cat. No.	Model	용량	Size (Ø×h, mm)	Neck id (mm)	PK Price ₩		CS Price ₩	
B08-46-621	T881080	125 ml	62×63	70	4ea	16,900	6pk	94,800
B08-46-623	T881090	250 ml	62×103		4ea	22,000	6pk	123,300
B08-46-625	T881100	500 ml	110×76	120	4ea	42,800	4pk	160,000
B08-46-627	T881111	1 L	110×135		4ea	56,400	4pk	210,600

- 재질 : PC Body, -135 ~ 135°C 내열, Autoclavable
- 투명 재질로 충격에 강함
- Cold Room이나 Refrigerator에서 사용 가능
- 요청에 의해 Blue, Green, Gray 제품도 공급 가능함
- Tarsons

