

**PP Container**  
**플라스틱 용기 (병 / 바이알)**

Cat. No.	Model	용량 (ml)	Size (몸통 od/캡 od×h)	눈금 (ml)	EA Price ₩	PK Price ₩
B08-660-001	4020-1806-13	30	31/37×44 mm	5 ~ 15	200	100ea 18,100
B08-660-005	4020-1906-13	50	44/50×42 mm	10 ~ 40	300	100ea 27,800
B08-660-009	4020-2006-13		34/40×70 mm		300	100ea 27,800
B08-660-013	4020-2106-13	80	44/51×61 mm	10 ~ 60	350	100ea 32,600
B08-660-017	4020-2206-13	150	60/53×81 mm	10~120	400	100ea 37,200
B08-660-021	4020-2306-13	250	59/65×100 mm	10~200	600	100ea 55,000

- 재질 : PP, -10 ~ 135°C 내열, Autoclavable
- Triple Molded Graduation
- LK Labkorea



**Sample Container, PP**  
**플라스틱 용기**

- 재질 : PP, -10 ~ 135°C 내열, Autoclavable
- Cap은 스크류 타입으로 체결이 간편함
- Rayhope®

**Low Form / 단형**

Cat. No.	Model	용량	Size	EA Price ₩
C21-126-064	LACE-01	100 ml	∅70×h50 mm	2,000
C21-126-070	LACE-02	200 ml	∅85×h60 mm	2,200
C21-126-076	LACE-03	400 ml	∅105×h75 mm	2,400
C21-126-082	LACE-04	900 ml	∅120×h85 mm	2,600



**High Form / 장형**

Cat. No.	Model	용량	Size	EA Price ₩
C21-126-088	LACE-05	250 ml	∅70×h95 mm	2,400
C21-126-094	LACE-06	450 ml	∅85×h110 mm	2,600
C21-126-100	LACE-07	900 ml	∅105×h150 mm	2,800
C21-126-106	LACE-08	1.35 L	∅120×h165 mm	2,900

**Plastic Container**  
**플라스틱 용기, 이너씰 및 캡 별매**

- 재질 : PP 300~1.3L , HDPE 1.5L 이상
- 1/4 회전 만에 개폐 가능
- Inner Seal 및 Cap 별도 구매
- Qorpak

**Container**

Cat. No.	Model	용량	재질	Neck Finish	Size (∅×h, mm)	EA Price ₩	PK Price ₩
B08-134-651	248735	300 ml	PP	96.8mm	99×55.4	5,700	90ea 475,800
B08-134-654	248736	500 ml			99×82.4	5,800	90ea 488,400
B08-134-657	248737	650 ml			99×104.4	6,200	90ea 516,700
B08-134-660	248738	1 L			99×157.4	6,800	90ea 570,200
B08-134-663	248740	1.3 L			99×197.4	7,800	60ea 432,800
B08-134-666	248742	1.5 L	HDPE	116mm	119×161	11,600	60ea 646,700
B08-134-669	248743	2 L			119×236	12,500	30ea 350,100
B08-134-672	248745	2.5 L			119×266	15,000	30ea 418,100



Container

**Inner Seal and Lid**

Cat. No.	Model	Style	사용 가능 용기	Material	EA Price ₩	PK Price ₩
B08-134-671	248746	Innner Seal	300 ml ~ 1.3 L	투명 PE	1,100	30ea 29,300
B08-134-674	248747		1.5 L ~ 2.5 L		2,300	30ea 62,000
B08-134-681	248748	Lid	300 ml ~ 1.3 L	Red PP	2,700	30ea 75,400
B08-134-683	248749		1.5 L ~ 2.5 L		5,600	30ea 154,600
B08-134-685	248750	Lid with Gasket	300 ml ~ 1.3L		3,200	30ea 89,400

- 248750은 Gasket 장착형



Inner Seal & Lid