

Sterile Storage Bottle, Square 사각 멸균 병

- 재질 : USP Class VI Class USP <661>, ISO10993
- Sterile : DNase-free, RNase-free, Non-pyrogenic
- 저온 및 자외선에 대한 내성이 높으며 파손 위험이 적음
- E-Beam Sterilization (전자빔 멸균)
- 버퍼, 배양액, pH에 민감한 샘플 보관용으로 알맞음

• LK Labkorea

PETG Bottle (for Pharmaceutical Industry)

Cat. No.	Model	용량	Size (Ø × h)	CS Price ₩	
B08-491-001	354111	30 ml	38 × 60.6 mm	40ea	113,500
B08-491-005	354511	60 ml	40 × 80.5 mm	48ea	124,800
B08-491-009	353511	125 ml	52.4 × 103 mm	48ea	141,800
B08-491-013	352511	250 ml	58.6 × 141 mm	48ea	239,600
B08-491-017	333511	500 ml	73 × 174 mm	24ea	151,000
B08-491-021	334511	1 L	92 × 215 mm	12ea	110,600

PET Bottle (for Laboratory)

Cat. No.	Model	용량	Size (Ø × h)	CS Price ₩	
B08-491-035	354611	60 ml	40 × 80.5 mm	48ea	73,700
B08-491-039	353611	125 ml	52.4 × 103 mm	48ea	83,300
B08-491-043	352611	250 ml	58.6 × 141 mm	48ea	127,600
B08-491-047	333001	500 ml	73 × 174 mm	24ea	100,700
B08-491-051	334001	1 L	92 × 215 mm	12ea	85,100



Sterile Media Bottle, Square PETG 사각 멸균 미디어 병

Cat. No.	Model	용량	Size (w × d × h)	Closure	Tray	PK Price ₩	
B08-46-301	T581410	30 ml	38.2 × 38.2 × 62.5 mm	20-415	24ea	96ea	118,400
B08-46-303	T581411	60 ml	40.4 × 40.4 × 81.1 mm	24-415		96ea	146,600
B08-46-305	T581412	125 ml	53.4 × 53.4 × 105 mm			48ea	103,700
B08-46-307	T581413	250 ml	60 × 60 × 143 mm	38-430	12ea	48ea	142,900
B08-46-309	T581414	500 ml	74.5 × 74.5 × 172 mm			24ea	118,400
B08-46-311	T581415	1 L	94 × 94 × 215 mm			24ea	196,500
B08-46-313	T581416	2 L	117 × 117 × 267 mm	53B	6ea	12ea	215,800

- 재질 : PETG, 75°C 내열
- USP Class VI, USP 37(661) & ISO10993-3을 준수함
- Leak-proof Seal & Gas Barrier
- Radiation Sterile : DNase / RNase Free, Non-cytotoxic, Non-pyrogenic, Animal Derived Component Free
- 뚜껑과 캡 사이에 Heat-shrink Band를 가지고 있음

• Tarsons



Sterile Diagnostic Bottle, Round PETG 멸균 원형 병, 5 ~ 20 ml

Cat. No.	Model	용량	Size (Ø × h)	Closure	Tray	PK Price ₩	
B08-46-490	T581701	5 ml	22.5 × 46.7 mm	20-415	20ea	100ea	75,300
B08-46-493	T581702	10 ml	24 × 57.5 mm			100ea	86,700
B08-46-496	T581703	20 ml	29.8 × 65.1 mm			100ea	116,700

- 재질 : PETG, 75°C 내열
- 기타 자세한 내용은 상단 Square PETG Bottle 참조
- 진단 테스트용 용기로 사용하기에 알맞음

• Tarsons



Sterile PET Sample Bottle PET 무균 검체병

Cat. No.	Model	용량	Size (Ø × h)	EA Price ₩	PK Price ₩	
B08-84-515	SSB0150	150 ml	55 × 88 mm	1,100	250ea	179,000
B08-84-525	SSB0250	250 ml	65 × 98 mm	1,300	150ea	130,000
B08-84-550	SSB0500	500 ml	75 × 134 mm	1,700	100ea	113,900
B08-84-585	SSB0850	850 ml	80(62) × h165 mm	2,700	50ea	89,600

- 재질 : PET Body, PE Lined PP Cap
- 개별 멸균 포장되어 있음

• LK Labkorea

