


I Bottle, PET 보급형 PET 병

- 뚜껑 재질 : PP (Polypropylene)
- Non-autoclavable
- Kyungshin Chem

Round PET Bottle

배송비 별도

Cat. No.	Model	용량	Size (Ø × h)	PK Price ₩	
B08-150-010	LKG102	25 ml	32 × 57 mm	650ea	195,800
B08-150-020	LKG110	100 ml	48 × 80 mm	252ea	132,300
B08-150-025	LKG120	200 ml	52 × 117 mm	245ea	147,600
B08-150-030	LKG125	250 ml	65 × 96 mm	216ea	150,600
B08-150-035	LKG130	300 ml	84 × 83 mm	128ea	118,900
B08-150-040	LKG150	500 ml	75 × 133 mm	100ea	79,200
B08-150-050	LKG850	850 ml	84 × 176 mm	84ea	74,500



Square PET Bottle

배송비 별도

Cat. No.	Model	용량	Size (w × d × h)	PK Price ₩	
B08-150-060	LKG210	120 ml	51 × 51 × 86 mm	210ea	119,300
B08-150-070	LKG220	200 ml	52 × 52 × 110 mm	245ea	147,600
B08-150-080	LKG235	350 ml	68 × 68 × 107 mm	150ea	116,100
B08-150-090	LKG250	500 ml	75 × 75 × 148 mm	120ea	95,100
B08-150-100	LKG295	950 ml	97 × 72 × 165 mm	100ea	104,100

I PET Bottle, Amber 보급형 갈색 PET 병

- 뚜껑 재질 : PP (Polypropylene)
- Non-autoclavable
- Kyungshin Chem

Amber Round PET Bottle

배송비 별도

Cat. No.	Model	용량	Size (Ø × h)	PK Price ₩	
B08-150-110	LKGA102	25 ml	32 × 57 mm	650ea	207,000
B08-150-120	LKGA110	100 ml	48 × 80 mm	252ea	136,600
B08-150-130	LKGA125	250 ml	65 × 96 mm	216ea	158,000
B08-150-140	LKGA150	500 ml	75 × 133 mm	100ea	83,500



Amber Square PET Bottle

배송비 별도

Cat. No.	Model	용량	Size (w × d × h)	PK Price ₩	
B08-150-160	LKGA210	120 ml	51 × 51 × 86 mm	210ea	124,700
B08-150-170	LKGA220	200 ml	52 × 52 × 110 mm	350ea	240,400
B08-150-180	LKGA235	350 ml	68 × 68 × 107 mm	150ea	129,000
B08-150-190	LKGA250	500 ml	75 × 75 × 148 mm	120ea	100,100
B08-150-200	LKGA295	950 ml	97 × 72 × 165 mm	100ea	109,300

I Water Collection Bottle, Sterilized PE 무균 채수병

Cat. No.	Model	용량	PK Price ₩	
B08-84-110	84WCB1N	1 L	150ea	110,000
B08-84-210	84WCB2N	2 L	100ea	94,900
B08-84-310	84WCB3N	4 L	65ea	93,400
B08-84-400	84WCB4	6 L	50ea	120,500

- Gamma-Ray 멸균 개별 포장
- Mold 눈금, Handle 부착형
- Flexible 재질로 운반에 편리
- 환경분야 샘플링에 적합
- LK Labkorea

