

## PFA Sample Bottle / Jar PFA 테프론 샘플 병 / 용기

Cat. No.	Model	용량	외부 Size (mm)	내부 Size (mm)	EA Price ₩
J01-460-031	A130-01	7 ml	∅ 22 x h 37	∅ 18.5 x h 32	109,300
J01-460-035	A130-05	60 ml	∅ 50 x h 47	∅ 45.5 x h 40	249,900
J01-460-039	A130-06	90 ml	∅ 50 x h 67	∅ 45.5 x h 58	295,800
J01-460-043	A130-07	120 ml	∅ 66 x h 53	∅ 60 x h 45	359,300
J01-460-047	A130-09	240 ml	∅ 66 x h 95	∅ 60 x h 86	489,900
J01-460-051	A130-11	500 ml	∅ 85 x h 120	∅ 80 x h 109	789,300
J01-460-055	A130-12	1 L	∅ 107 x h 151	∅ 100 x h 139	1,521,000

- 재질 : PTFE , Micro Porous PTFE , -200 °C to 250°C
- 투명한 테프론 재질로 최상의 내화학성
- BOLA



## PFA Jar with Ground Joint PFA 테프론 자, 상부 조인트 ST 29/32

Cat. No.	Model	용량	Size	i.d	EA Price ₩
J01-460-081	A159-24	1L	∅ 107 x h193mm	∅ 100mm	1,869,400

- 재질 : PFA Body, PFA Cap with PTFE Joint
- 상부 Joint : ST 29/32
- 투명한 테프론 재질로 최상의 내화학성
- 용도 : Evaporation Vessel, Reaction Vessel
- BOLA



## PFA Jar with Tubing Port PFA 테프론 자, 튜빙 연결 피팅

Cat. No.	Model	용량	외부 Size	i.d	EA Price ₩
J01-460-011	A131-12	120 ml	∅ 66 x h 74mm	∅ 60mm	605,900
J01-460-015	A131-14	240 ml	∅ 66 x h 116mm	∅ 60mm	736,600

- 재질 : PTFE , -200 °C to 250°C 내열
- 2개의 Tubing 연결 피팅 : 1/4"(0.64mm) Tubing 용
- 투명한 테프론 재질로 최상의 내화학성
- BOLA



## Aluminum Bottle with PP Screw Cap 순수 알루미늄 바틀

Cat. No.	Model	용량	Size	EA Price ₩
B08-100-011	9072168	38 ml	∅35×h70 mm	4,100
B08-100-015	9072169	60 ml	∅40×h88 mm	4,500
B08-100-019	9072170	120 ml	∅57×h97 mm	5,400
B08-100-023	9072171	300 ml	∅57×h160 mm	7,200
B08-100-027	9072172	600 ml	∅74×h192 mm	10,500
B08-100-031	9072173	1200 ml	∅88×h248 mm	12,400

- 재질 : Pure Aluminium (AL 99.5 Material)
- PP Screw Cap with Internal Aluminium Seal
- 변형되지 않은 상태로 보관해야 하는 샘플을 저장, 운반하는데 적합
- Wenk Labtec

