

Stainless Steel Bottle 스테인레스 바틀

Cat. No.	Model	용량	Size	EA Price ₩
B08-695-101	92-0041	200 ml	Ø60×h124 mm	310,600
B08-695-103	92-0042	500 ml	Ø80×h165 mm	332,700
B08-695-105	92-0043	1 L	Ø100×h206 mm	353,000
B08-695-107	92-0044	2 L	Ø120×h258 mm	385,300
B08-695-109	92-0045	5 L	Ø180×h319 mm	492,200
B08-695-111	92-0046	10 L	Ø220×h386 mm	599,200

- 재질 : Stainless Steel 304
- Thickness : 1.0 mm
- Cap 재질 : SUS303 (내부는 PTFE Packing)
- Sansyo



Wide Mouth Stainless Steel Bottle 스테인레스 광구병

Cat. No.	Model	용량	Size	EA Price ₩
B08-695-113	92-0031	2 L	Ø120×h218 mm	1,011,600
B08-695-115	92-0032	5 L	Ø160×h311 mm	1,098,200
B08-695-117	92-0033	10 L	Ø210×h376 mm	1,183,000

- 재질 : Stainless Steel 304, PTFE 패킹
- 입구 Size : Ø91 mm
- Sansyo



Stainless Steel Bottle, Clamp Type 클램프형 스테인레스 바틀

Cat. No.	Model	용량	Size	EA Price ₩
B08-695-119	92-0034	2 L	Ø120×h239 mm	779,000
B08-695-121	92-0035	5 L	Ø160×h324 mm	875,800
B08-695-123	92-0036	10 L	Ø210×h374 mm	972,500

- 재질 : Stainless Steel 304, Silicon Rubber 패킹
- 직경이 Ø97.6 mm 넓기 때문에 세척이 용이
- 클램프로 조이는 방식으로 새기 쉬운 내용물의 저장 및 운반에 적합
- Sansyo



Straight Side Stainless Steel Bottle 스테인레스 바틀

Cat. No.	Model	용량	Size	EA Price ₩
B08-695-125	92-1895	100 ml	Ø47.8×h62 mm	578,300
B08-695-127	92-1896	200 ml	Ø59.5×h85 mm	683,700
B08-695-129	92-1897	500 ml	Ø72.3×h135 mm	789,300
B08-695-131	92-1898	1 L	Ø97.6×h150 mm	980,200

- 재질 : Stainless Steel 304 (Body), Stainless Steel 316L (Cap), Silicon Rubber 패킹
- Cap에 Handle이 달려있어 개폐가 용이
- Sansyo

