

I Dropping Bottle, Screw Cap 스포이드 병, 스크류 캡

- Sodalime Glass Bottle and Black Screwcap / Dropper
- Autoclavable • LK Labkorea

Round-type, PP Cap, 10 ~ 30ml

Cat. No.	Model	용량	Height	Description	EA Price ₩
B08-53-260	L120200	10 ml	75 mm	Clear	1,300

Square-type, Phenolic Cap, 30 ~ 100ml

Cat. No.	Model	용량	Height	Description	EA Price ₩
B08-53-280	L120105	30 ml	104 mm	Clear	1,000
B08-53-284	L120106	50 ml	114 mm		1,100
B08-53-288	L120107	100 ml	140 mm		1,400
B08-53-296	L120109	50 ml	114 mm	Amber	1,200
B08-53-300	L120110	100 ml	140 mm		1,400



I Dropping Bottle 드라핑 바틀

Cat. No.	Model	용량	Neck	Size (Ø×h)	EA Price ₩
B08-59-411	125216	20 ml	GL14	31×88 mm	1,200
B08-59-413	125224	30 ml		34×96 mm	1,300
B08-59-415	125228	50 ml	GL18	39×115 mm	1,600
B08-59-417	125238	100 ml		45×136 mm	1,700
B08-59-419	125248	250 ml	GL25	59×170 mm	2,700
B08-59-421	125254	500 ml		75×209 mm	4,200
B08-59-423	125262	1 L	GL28	94×240 mm	9,800

- 재질 : LDPE, -100 ~ 80°C 내열, Non-autoclavable
- Dropping Nozzle과 Spout Cap 포함
- Brand GmbH



I PTFE Dropping Bottle PTFE 테프론 드라핑 바틀

Cat. No.	Model	용량	Size	EA Price ₩
B08-450-190	014.025	25 ml	Ø33×h80 mm	72,000
B08-450-195	014.050	50 ml	Ø43×h100 mm	78,300

- 재질 : PTFE, -200 ~ 280°C 내열
- Dropping Nozzle, PTFE Cap 포함
- Flexible
- Cowie



I Dropping Bottle / Dispensing Bottle 드라핑 바틀

Cat. No.	Model	용량	Cap	PK Price ₩	CS Price ₩
B08-46-321	T620030	15 ml	20 mm	12ea 10,100	6pk 56,500
B08-46-323	T620040	30 ml		12ea 12,200	6pk 68,300
B08-46-325	T620050	60 ml	24 mm	12ea 15,000	4pk 55,900
B08-46-327	T620060	125 ml		12ea 19,500	4pk 72,700
B08-46-329	T620070	250 ml		6ea 14,600	6pk 81,900

- 재질 : LDPE Body, PP Closure
- USP Class VI에 따른 Medical Grade Resin 사용
- Tarsons

