

Lab Spray Bottle / Automizer, HDPE

분무기

Cat. No.	Model	용량	EA Price ₩	PK Price ₩
B08-59-881	144180	400 ml	18,500	5ea 86,400
B08-59-883	144190	850 ml	19,900	5ea 92,900

- 재질 : HDPE
- Nozzle를 돌려 Spray 되어지는 정도를 조절 할 수 있음
- Brand GmBH



Trigger Spray Bottle, HDPE & Nylon

분무기

Cat. No.	Model	용량	Bottle	EA Price ₩
B08-134-107	PLC03499	650 ml	HDPE / Natural	12,300
B08-134-108	PLC04489	960 ml	Nylon / Natural	10,600
B08-134-109	PLC05489		Nylon / Amber	11,000

- 재질 : HDPE, -100 ~ 120°C 내열, Non-autoclavable
- HDPE From Lined Cap • 내화학성 분무기 • Qorpak



Wash Bottle with Adjustable Sprayer

분무기

Cat. No.	Model	용량	Description	EA Price ₩
B08-147-050	3345 ①	240 ml	Adjustable Spray	13,500
B08-147-100	3340 ②	1 L		65,100

- PE(Polyethylene) 내화학성 분무기
- Side에 또다른 Cap이 있어 용액을 채우기에 용이
- Cole-Parmer



Autoclavable Spray Bottle

멸균용 분무기

Cat. No.	Model	용량	용도	EA Price ₩
B08-215-151	885	500 ml	Ethanol, Iodine based Disinfectant 외	40,400
B08-215-153	887-03IP		Isopropyl alcohol (IPA), Ethanol 외	48,400

- 재질 : PP Bottle & Nozzle, Silicon Hose • 멸균 가능 : 120°C, 20min
- For Pharmaceutical Production, QC Testing and Clinical Laboratory
- Flexible Hose with Weight with a Filter
- Universal



PE Pump Vaporizer

PE 스프레이형 분무기

Cat. No.	Model	용량	분주량	Size	EA Price ₩
B08-59-888	LSB020	20ml	0.12ml	∅29.5×h62mm	1,000
B08-59-891	LSB050	50ml		∅32.2×h87mm	1,200
B08-59-893	LSB100	100ml		∅37.5×h118mm	1,300
B08-59-896	LSB250	250ml		∅49×h164mm	1,500

- 재질 : LDPE (20~50ml), HDPE (100 ~250ml) Body / PP Closure
- 스프레이식 분사형 • Bottle 및 Closure 반투명
- LabPlasti®



Pumping Dispenser Bottle, HDPE

HDPE 분주 바들

Cat. No.	Model	용량	Size	EA Price ₩
B08-20-401	LDB150	150ml	∅42.5×h133mm	1,400
B08-20-405	LDB250	250ml	∅49×h164mm	1,500
B08-20-409	LDB500	500ml	∅63×h200mm	1,700
B08-20-413	LDB1000	1 L	∅80×h247mm	2,400

- 재질 : HDPE Body, PP Closure
- Gel과 같이 점도가 있는 용액을 담고 분주할 때 쓰임 • LabPlasti®

