

I Narrow Mouth Handled Jug Bottle, HDPE

세구 핸들 대용량 병

Cat. No.	Model	용량	Size	EA Price ₩
B08-134-001	PLC03504	1.9 L	Ø121×h254 mm	8,800
B08-134-003	PLC03518	3.8 L	Ø157×h298 mm	8,000

- 재질 : HDPE, PP Cap with PE Form Liner
- Cap Size : 38-400
- Qorpak



I Wide Mouth Handled Jug Bottle, HDPE

광구 핸들 대용량 병

Cat. No.	Model	용량	Size	EA Price ₩
B08-134-007	PLC06548	3.8 L	Ø159×h232 mm	11,800

- 재질 : HDPE, PP Cap with PE Form Liner
- Cap Size : 89-400
- Qorpak



I Narrow Mouth Large Capacity Bottle

대용량 PE 세구병

Cat. No.	Model	용량	Size (Neck id)	EA Price ₩
B08-59-711	130870	5 L	Ø165×h335 mm (45 mm)	44,600
B08-59-713	130874	10 L	Ø210×h420 mm (54 mm)	63,300

- 재질 : HDPE, -100 ~ 120°C 내열, Non-autoclavable
- Carrying Handle를 가지고 있어 이동이 편리함
- Brand GmBH



I Wide Mouth Large Capacity Bottle

대용량 PE 광구병

Cat. No.	Model	용량	Size (Neck id)	EA Price ₩
B08-59-731	130970	5 L	Ø165×h320 mm (85 mm)	44,800
B08-59-735	130974	10 L	Ø210×h395 mm (120 mm)	63,700

- 재질 : HDPE, -100 ~ 120°C 내열, Non-autoclavable
- Carrying Handle를 가지고 있어 이동이 편리함
- Brand GmBH

