

Carboy with Stopcock, PP PP 증류수 통 [병]

Cat. No.	Model	용량	Size	EA Price ₩
B08-46-413	T583210	5 L	Ø168×h357 mm	52,000
B08-46-420	T583280	10 L	Ø249×h390 mm	80,200
B08-46-427	T583290	20 L	Ø285×h535 mm	102,400
B08-46-430	T583300	50 L	Ø376×h682 mm	179,000

- 재질 : PP, -10 ~ 135°C 내열, Autoclavable
- USP Class VI에 따른 Medical Grade Resin 사용
- Leakproof Closure, Graduation
- Spigot Closure 기본 포함
- Cap Size : 83 mm (id65 mm), 단, 5 L는 53 mm (Neck id : 43.5 mm)임
- Tarsons

Optional Accessory

Cat. No.	Model	Description	EA Price ₩
B08-46-253	T583217	Filing / Venting Closure for 83 mm Cap	62,500
B08-46-371	T583212	Spare Cap, 83 mm	9,900
B08-46-406	T583211	Spare Cock / Spigot	27,700



Filing / Venting Closure

Quick Filling & Venting Closure for Carboy 카보이 퀵 컨넥터 뚜껑, 83B

Cat. No.	Model	Cap	Fitting	EA Price ₩
B08-76-643	L215821	83B	1/4" × 2	264,700
B08-76-646	L215823		3/8" × 2	268,100
B08-76-649	L215831		1/4" × 3	360,000
B08-76-652	L215833		3/8" × 3	375,800

- 재질 : PP Cap
- Quick Disconnect Fitting
- Inside & Outside Barbed Tubulation
- LK Labkorea



Air Venting Filter는 419 Page 참조



Carboy with Stopcock, Rectangular, PP PP 사각 증류수 통 [병]

Cat. No.	Model	용량	Size	EA Price ₩
B08-46-434	T683240	5 L	w227×d144×h342 mm	60,800
B08-46-441	T683250	10 L	w284×d176×h400 mm	83,700
B08-46-448	T683260	20 L	w342×d207×h482 mm	117,600
B08-46-406	T583211	Spare Cock / Spigot		27,700

- 재질 : PP, -10 ~ 135°C 내열, Autoclavable
- USP Class VI에 따른 Medical Grade Resin 사용
- 공간 절약형의 사각 증류수 통
- Cap Size : 100 mm (Neck Id : 85.5 mm)
- Leakproof Closure, Graduation
- Spigot Closure 기본 포함, Stainless Steel Handle
- Tarsons

