


### Gas Washing Bottle 스크류 캡 가스 세척병, LukeGL®

Cat. No.	Model	용량	Description	EA Price ₩
B08-83-168	530038	100 ml	With Filter Disc P2, 40 ~ 100 um ①	85,700
B08-83-176	530048	250 ml		92,300
B08-83-184	530058	500 ml		99,200
B08-83-192	530138	100 ml	Without Filter Disc ②	69,900
B08-83-200	530148	250 ml		72,300
B08-83-208	530158	500 ml		75,900

- Borosilicate Glass 3.3
- with Screw Connection GL14
- Porosity 규격 선택 가능
  - P0 (160 ~ 250 um), P1 (100 ~ 160 um), P2 (40 ~ 100 um), P3 (16 ~ 40 um), P4 (10 ~ 16 um), P5 (1 ~ 1.6 um)
- LukeGL®



### Gas Washing Bottle 가스 세척병, LukeGL®

- Borosilicate Glass 3.3
- Bottle and Head to DIN 12596
- Without Filter Disc
- LukeGL®

#### Hose Nozzle 0~ Ø11mm

Cat. No.	Model	용량	Conical Joint	EA Price ₩
B08-83-008	5500537	100 ml	29/32	43,400
B08-83-016	5500549	250 ml		52,700
B08-83-024	5500558	500 ml	34/35	66,300
B08-83-032	5500570	1 L	45/40	89,500

#### Thread GL14

Cat. No.	Model	용량	Conical Joint	EA Price ₩
B08-83-040	55005G1	100 ml	29/32	55,400
B08-83-048	55005G2	250 ml		63,900
B08-83-056	55005G3	500 ml	34/35	78,300
B08-83-064	55005G4	1 L	45/40	101,600



### Gas Washing Bottle 필터 부착형 가스 세척병, LukeGL®

Cat. No.	Model	용량	Porosity	Conical Joint	EA Price ₩
B08-83-072	5510037	100 ml	160 ~ 250 um P0	29/32	55,500
B08-83-080	5510049	250 ml		34/35	62,700
B08-83-088	5510058	500 ml		45/40	86,400
B08-83-096	5510070	1 L		45/40	101,100
B08-83-104	5510137	100 ml	100 ~ 160 um P1	29/32	55,500
B08-83-112	5510149	250 ml		34/35	62,700
B08-83-120	5510158	500 ml		45/40	86,400
B08-83-128	5510170	1 L		45/40	101,100
B08-83-136	5510237	100 ml	40 ~ 100 um P2	29/32	55,500
B08-83-144	5510249	250 ml		34/35	62,700
B08-83-152	5510258	500 ml		45/40	86,400
B08-83-160	5510270	1 L		45/40	101,100

- Borosilicate Glass 3.3
- Large Disc Provided Greater 용량
- Hose Nozzle Ø0 ~ 11 mm
- 그외 Porosity 규격 선택가능
  - P3 (16 ~ 40 um), P4 (10 ~ 16 um), P5 (1 ~ 1.6 um)
- LukeGL®

