

Density Bottle w/ Thermometer / Pycnometer 비중병, 온도계 부착형

Cat. No.	Model	용량	Tolerance	EA Price ₩
B08-83-235	5640013	10 ml	±0.001 ml	338,300
B08-83-243	5640022	25 ml		346,300
B08-83-251	5640028	50ml		357,600
B08-83-259	5640037	100ml		376,300

- Borosilicate Glass 3.3
- 조인트 온도계 부착형으로서 ST10/19 with Side Tube ST7/16 Cap
- Thermometer : 5 ~ 35°C (0.2°C 눈금)
- Calibrated to Contain (TC, In)
- 용도 : 분말시료의 비중을 습식방법으로 측정하는 비중병
- Lenz



Pycnometer with Thermometer 온도계 부착형 비중병, Kimble®

Cat. No.	Model	용량	Height	EA Price ₩
B08-137-296	15123N-10	10 ml	168 mm	407,200
B08-137-300	15123N-25	25 ml	175 mm	407,200
B08-137-304	15123N-50	50 ml	186 mm	407,200

- Kimax® Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- 액체 또는 고체의 비중 측정 및 비교용
- Non-Mercury Filled Thermometer 부착
- 조인트형 온도계 : 14 ~ 38°C (0.2°C 눈금)
- 5/12 Side Tube Cap 포함
- DWK Lifescience / Kimble®



Density Bottle / Pycnometer 비중병, Guy-Lussac, 보증서 포함

Cat. No.	Model	용량	Tolerance	EA Price ₩
B08-52-241	43305	5 ml	± 1.00 ml	93,700
B08-52-243	43308	10 ml		96,900
B08-52-247	43320	25 ml	± 2.00 ml	106,100
B08-52-230	43328	50 ml	± 3.00 ml	121,100
B08-52-233	43338	100 ml		165,200

- Borosilicate Glass 3.3, USP Type I, Class A
- Complies with ISO 3507, Guy-Lussac type
- Stopper : NS 10/19 with Capillary
- Calibrated to Contain (TC, In)
- 보증서 (Batch Certificate) 포함
- Brand GmBH



Specific Gravity Bottle / Pycnometer 비중병

Cat. No.	Model	용량	Joint	EA Price ₩
B08-83-216	5660006	500 ml	45/27	151,000
B08-83-224	5660007	1 L	60/37	170,600

- Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- For Density Determination
- Consists of Wide-neck Bottle, Conical Shoulder
- With NS Joint and Pycnometer Head
- Lenz

