

**Wide Mouth Wash Bottle, Side Delivery Tube**  
**광구 일체형 세척병, LDPE**

Cat. No.	Model	용량	Cap	EA Price ₩
B08-46-901	T561100	125 ml	24 mm	1,900
B08-46-903	T561110	250 ml	38 mm	3,300
B08-46-905	T561120	500 ml		3,900
B08-46-907	T561130	750 ml		4,700

- 재질 : LDPE, -100 ~ 80°C 내열
- 완전 일체형 Washing Bottle
- Tarsons



**Safety Vented Wash Bottle, Side Delivery Tube**  
**벤트 일체형 안전 세척병, PTFE 벤트 캡 포함**

Cat. No.	Model	용량	Cap 색상	EA Price ₩
B08-46-881	T561120VN	500 ml	Natural	6,300
B08-46-883	T561120VR		Red	6,300
B08-46-885	T561120VY		Yellow	6,300
B08-46-887	T561120VG		Green	6,300
B08-46-889	T561120VB		Blue	6,300
B08-46-891	T561130VN	750 ml	Natural	8,100
B08-46-893	T561130VR		Red	8,100
B08-46-895	T561130VY		Yellow	8,100
B08-46-897	T561130VG		Green	8,100
B08-46-899	T561130VB		Blue	8,100

- 재질 : LDPE, PP Closure
- 벤트용 완전 일체형 세척병
- 마개 상단에 Solvent Venting용 PTFE Membrane이 부착되어 있음
- Tarsons



**PFA Washing Bottle**  
**PFA 테프론 세척병**

- 재질 : PFA, -196 ~ 250°C 내열
- Flon Industry

**Wide-Mouth, 광구 세척병 ①**

Cat. No.	Model	용량	EA Price ₩
B08-64-922	F1011-01	100 ml	118,800
B08-64-925	F1011-02	250 ml	161,200
B08-64-928	F1011-03	500 ml	254,600
B08-64-931	F1011-04	1 L	446,400

**Narrow-Mouth, 세구 세척병 ②**

Cat. No.	Model	용량	EA Price ₩
B08-64-962	F1010-01	100 ml	98,400
B08-64-965	F1010-02	250 ml	129,000
B08-64-968	F1010-03	500 ml	176,500
B08-64-971	F1010-04	1 L	404,000

**PFA Wash Bottle Cap Only, 세척병용 캡 ③**

Cat. No.	Model	용량	EA Price ₩
B08-64-981	F1014-001	100 ml 용	34,000
B08-64-983	F1014-002	250 ml 용	34,000
B08-64-985	F1014-003	500 ml 용	34,000
B08-64-987	F1014-004	1 L 용	34,000

