



Chemical Hazard Label 케미컬 라벨 테이프

Cat. No.	Model	Size (w×d)	규격	EA Price ₩
B08-218-053	SBH-30	63.5×101.6 mm	1000/Roll	239,800
B08-218-055	SBH-31	50.8×50.8 mm		107,200
B08-218-057	SBH-32	25.4×25.4 mm		64,200
B08-218-101	UPCR-7002	25.4×76.2 mm	500/Roll	79,100
B08-218-105	UPCR-7003	25.4×76.2 mm		79,100
B08-218-131	UPCR-7036	50.8×76.2 mm		93,000


• Shamrock




UPCR-7036




SBH-30




UPCR-7002



SBH-31



UPCR-7003



SBH-32

Safety Solvent Bottle Carrier 안전 솔벤트 바틀 캐리어, Tote Style

Cat. No.	Model	Type	색상	EA Price ₩
B08-134-970	235349	500 ml & 1 L 용	Red	135,200
B08-134-974	235350		Black	135,200
B08-134-978	235345	2.5 L & 4 L 용	Red	216,900
B08-134-982	235346		White	216,900
B08-134-986	235348		Black	216,900

- 재질 : Thermoplastic Elastomer, -45 ~ 110°C 내열, Autoclavable
- Most Aqueous Solutions, Acids & Bases에 내화학성은 가짐
- Bacteria & Fungi Growth에 내성이 있음
- 손잡이 일체형으로 본체와 하나의 Mold로 제작
- Excellent Tensile Strength
- Liquid 와 Ice 의 이송용으로도 사용 가능
- 캐리어 하부에 Centering Ring이 있어 Bottle과 벽과의 거리를 유지 시킴
- Recessed Marking Area at Top Outer Side of Carrier
- 500 ml & 1 L 용 Size : Ø114 x h165 mm (Bucket), 전체 높이 254 mm
- 2.5 L & 4 L 용 Size : Ø190 x h266 mm (Bucket), 전체 높이 406 mm
- Qorpak



Safety Solvent Bottle Carrier 안전 솔벤트 바틀 캐리어, Cover Style

- 재질 : Epoxy Coated PE, Heavy-gauge Steel Wire Bail
- Wide Ribs in Both the Sides and on The Bottom to Support and Cushion Bottles
- Cone Style Lid : 16mm Opening in the Top for Insertion of Tubing
- Simplifies the Transfer of Liquids without Pouring
- Truncated Style Lid : Pouring without the Removal of the Lid
- Bel-Art

Cone Style, 16 mm Opening for Tubing Inserting ①

Cat. No.	Model	적용 Bottle	Size	EA Price ₩
B08-11-430	F16957-0000	2.5 L	Ø180×h350 mm	260,900
B08-11-434	F16958-0000	4 L	Ø200×h370 mm	309,500

Truncated Style, for Pouring ②

Cat. No.	Model	적용 Bottle	Size	EA Price ₩
B08-11-438	F16959-0000	2.5 L	Ø180×h320 mm	260,900
B08-11-442	F16960-0000	4 L	Ø200×h350 mm	307,800

