

Air Tight Box, PE 플라스틱 박스

Cat. No.	Model	용량	Size (w×d×h)	EA Price ₩
B09-20-325	ATB0500	500 ml	103×103×64 mm	6,100
B09-20-329	ATB1000	1 L	103×103×127 mm	9,500
B09-20-333	ATB1001	1 L	208×103×64 mm	11,700
B09-20-337	ATB1500	1.5 L	208×103×94 mm	15,700
B09-20-341	ATB2000	2 L	208×208×64 mm	21,300
B09-20-345	ATB3500	3.3 L	208×208×94 mm	22,300

- 재질 : PE, -40 ~ 95°C 내열
- 샘플의 저장, 이송 등에 사용되는 다용도 사각 박스
- LabPlasti®



PP Rectangular Air Tight Box PP 4각 밀폐용기

Cat. No.	Model	용량	Size (w×d×h)	EA Price ₩
B09-735-005	HPL805	180 ml	108×88×48 mm	2,700
B09-735-010	HPL806	350 ml	135×102×52 mm	3,300
B09-735-015	HPL807	470 ml	135×102×68 mm	4,000
B09-735-020	HPL815	550 ml	181×128×52 mm	4,900
B09-735-040	HPL808	850 ml	135×102×118 mm	4,500
B09-735-055	HPL817	1 L	205×134×69 mm	6,300
B09-735-095	HPL819	2 L	135×102×282 mm	8,100
B09-735-120	HPL827P	4.5 L	248×180×150 mm	22,200

- 재질 : PP, -10 ~ 135°C 내열, Autoclavable
- HPL827P 제품은 뚜껑 갈색임
- Lock & Lock



Clear Storage Box, PC 다용도 플라스틱 투명 박스

Cat. No.	Model	용량	Size (w×d×h)	EA Price ₩
B09-517-007	PCB007	7 L	300×450×90 mm	39,100
B09-517-012	PCB012	12 L	300×450×150 mm	40,000
B09-517-020	PCB020	20 L	300×450×230 mm	43,400
B09-517-030	PCB030	30 L	450×660×150 mm	60,800
B09-517-050	PCB050	50 L	450×660×230 mm	71,600
B09-517-065	PCB065	65 L	450×660×300 mm	93,300
B09-517-080	PCB080	80 L	450×660×380 mm	119,300

- 재질 : PC, -30 ~ 108°C 내열
- LK Labkorea



Cryo Storage Box, LDPE 냉동 보관 박스

Cat. No.	Model	Capacity	Size (mm)	EA Price ₩
B09-735-301	309040	120 ml	w100×d67×h38	3,900
B09-735-311	309110	180 ml	w100×d67×h55	5,200
B09-735-321	309280	300 ml	w138×d100×h38	4,700
B09-735-331	309350	450 ml	w138×d100×h55	6,400
B09-735-341	309420	650 ml	w206×d138×h38	6,900
B09-735-351	309590	1 L	w206×d138×h55	9,000
B09-735-361	324960	1.3 L	w285×d105×h72	9,900
B09-735-371	324890	3 L	w288×d210×h72	11,800

- 재질 : Body - SAN, Cover - LDPE
- Powder, Filter, Chromatography Sheet, Stopper, 냉동 보관 및 비휘발성 액체를 단기 보관 하기에 적합
- LK Labkorea

