

PTFE Cock Buret

테프론 콕 뷰렛, Class AS + Batch 보증서

- Borosilicate Glass 3.3
- Complied with DIN EN ISO 385
- 보증서 (Batch Certificate) 포함
- Techno Glass

Clear Buret / 투명 뷰렛

Cat. No.	Model	용량	눈금 (ml)	Tol. (ml)	EA Price ₩
B13-74-523	1567/AS-TK-10	10 ml	0.02	±0.02	44,500
B13-74-530	1567/AS-TK-25	25 ml	0.05	±0.03	42,500
B13-74-537	1567/AS-TK-50	50 ml	0.10	±0.05	43,400
B13-74-544	1567/AS-TK-100	100 ml	0.20	±0.10	56,100

Amber Buret / 갈색 뷰렛

Cat. No.	Model	용량	눈금 (ml)	Tol. (ml)	EA Price ₩
B13-74-556	1567/AS-TK-10A	10 ml	0.02	±0.02	71,800
B13-74-564	1567/AS-TK-25A	25 ml	0.05	±0.03	69,900
B13-74-572	1567/AS-TK-50A	50 ml	0.10	±0.05	70,700
B13-74-576	1567/AS-TK-100A	100 ml	0.20	±0.10	83,500

PTFE Cock Buret with Funnel Top

펀넬 탑 테프론 콕 뷰렛, Class AS + Batch 보증서

Cat. No.	Model	용량	눈금 (ml)	Tol. (ml)	EA Price ₩
B13-74-623	1567/AS-TKK-10	10 ml	0.02	±0.02	45,300
B13-74-630	1567/AS-TKK-25	25 ml	0.05	±0.03	43,100
B13-74-637	1567/AS-TKK-50	50 ml	0.10	±0.05	44,200
B13-74-644	1567/AS-TKK-100	100 ml	0.20	±0.10	56,900

- Borosilicate Glass 3.3
- Complied with DIN EN ISO 385
- 보증서 (Batch Certificate) 포함
- Techno Glass

Compact Buret, Modular

분리형 뷰렛

- 구성 : Glass Buret Bulb & Tip, PTFE Buret Stopcock
- Complied with DIN EN ISO 385, Calibrated to deliver (TD, Ex)
- Module형으로 PTFE Cock 및 Tip을 교체 할 수 있음
- Brand GmBH

Clear Buret, Class AS / 투명 분리형 뷰렛, Batch 보증서

Cat. No.	Model	용량	눈금 (ml)	Tol. (ml)	EA Price ₩
B13-59-721	13913	10 ml	0.02	±0.02	128,200
B13-59-725	13916	25 ml	0.05	±0.03	128,200
B13-59-729	13918	50 ml	0.1	±0.05	128,200

- Borosilicate Glass 3.3
- 보증서 (Batch Certificate) 포함
- BLAUBRAND®

Amber Buret, Class B / 갈색 분리형 뷰렛

Cat. No.	Model	용량	눈금 (ml)	Tol. (ml)	EA Price ₩
B13-59-735	13926	25 ml	0.1	±0.08	150,500
B13-59-739	13928	50 ml		±0.10	150,500

- Borosilicate Glass 5.0
- SILBERBRAND®

Spare PTFE Buret Stopcock

Cat. No.	Model	Description	EA Price ₩
B13-59-761	11805	PTFE Cock only, without Tip	70,100

- 모든 Compact Buret과 같이 사용

Spare Glass Buret Tip

Cat. No.	Model	색상	Type	EA Price ₩
B13-59-772	11510	Clear	10 ml Buret 용	11,800
B13-59-774	11511		25 ml Buret 용	11,800
B13-59-776	11512		50 ml Buret 용	11,800
B13-59-778	11515	Amber	25 ~ 50 ml Buret 용	14,000

- Borosilicate Glass 3.3

