


Serialized & Certified Buret, 3-way Stopcock ASTM 테프론 콕크 뷰렛, Class A + 개별 보증서

Cat. No.	Model	용량	눈금 (ml)	Tol. (ml)	EA Price ₩
B13-137-111	17115F-5	5 ml	0.01	±0.01	691,100
B13-137-115	17115F-10	10 ml	0.02	±0.02	613,000
B13-137-119	17052F-25	25 ml	0.1	±0.03	352,800
B13-137-123	17052F-50	50 ml	0.1	±0.05	371,700
B13-137-127	17052F-100	100 ml	0.2	±0.1	476,400

- Kimax® Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- 개별 시리얼 넘버가 Marking 되어 있고 USP 개별 보증서 포함
- Filling Tube 는 id1/4" 튜브에 맞맞음
- Bore 2 mm, 3-way Stopcock
- ASTM E287 Class A, ASTM E694, E542 - USP Grade
- DWK Lifescience / Kimble®



PTFE Needle Valve Buret 니들밸브형 투명 뷰렛, Class AS + Batch 보증서

Cat. No.	Model	용량	눈금 (ml)	Tol. (ml)	EA Price ₩
B13-59-521	12484	10 ml	0.02	±0.02	83,100
B13-59-525	12487	25 ml	0.05	±0.03	83,100
B13-59-529	12488	50 ml	0.1	±0.05	83,100

- Borosilicate Glass 3.3
- Complied with DIN EN ISO 385, Calibrated to deliver (TD, Ex)
- 보증서 (Batch Certificate) 포함
- 길이 : 약 800 mm, Schellbach Strips
- Brand GmBH



Glass Cock Buret 유리 콕크 뷰렛, Class AS + Batch 보증서

Cat. No.	Model	용량	눈금 (ml)	Tol. (ml)	EA Price ₩
B13-74-601	1567/AS-K-10	10 ml	0.02	±0.02	38,200
B13-74-603	1567/AS-K-25	25 ml	0.05	±0.03	36,300
B13-74-605	1567/AS-K-50	50 ml	0.10	±0.05	37,100
B13-74-607	1567/AS-K-100	100 ml	0.20	±0.10	49,900

- Borosilicate Glass 3.3
- Complied with DIN EN ISO 385
- 보증서 (Batch Certificate) 포함
- Techno Glass



Plastic Buret, PMMA PMMA 플라스틱 뷰렛

Cat. No.	Model	용량	EA Price ₩
B13-503-001	10330301	25 ml	13,300
B13-503-011	10330302	50 ml	14,800
B13-503-015	10330303	100 ml	24,000

- 재질 : Polymethyl Methacrylate (PMMA) / TPX, PTFE Cock
- 사용조건 : up to 5% wt/wt of Acids or 30% Base
- pH Range : 1 ~ 14
- LabPlasti®

