

## Micro Buret with Needle Valve 니들 밸브형 마이크로 뷰렛, Class AS + Batch 보증서

- Borosilicate Glass 3.3
- Complied with DIN EN ISO 385, Calibrated to deliver (TD, Ex)
- 보증서(Batch Certificate) 포함
- PTFE Stopcock인 중간밸브를 가지고 있음
- Brand GmBH

### Straight Needle Valve

Cat. No.	Model	Capacity (ml)	눈금 (ml)	Tol. (ml)	Overall Length(mm)	EA Price ₩
B13-59-400	24265	2	0.01	±0.01	600	313,400
B13-59-405	24267	5				318,900
B13-59-410	24269	10	0.02	±0.02	820	318,900

### Lateral Needle Valve

Cat. No.	Model	Capacity (ml)	눈금 (ml)	Tol. (ml)	Overall Length(mm)	EA Price ₩
B13-59-415	24595	2	0.01	±0.01	660	399,300
B13-59-425	24597	5				441,100
B13-59-420	24599	10	0.02	±0.02	900	441,100


- 하부 받침대 포함



## Micro Buret with PTFE Cock 마이크로 뷰렛, Class AS + Batch 보증서

Cat. No.	Model	용량	눈금	EA Price ₩
B13-74-125	436063314	1 ml	0.01 ml	114,900
B13-74-129	436063316	2 ml		115,200
B13-74-133	436063517	3 ml	0.02 ml	115,400
B13-74-137	436063817	5 ml		117,100

- Borosilicate Glass 3.3
- 시액 주입기 및 2개의 PTFE Stopcock
- 보증서 (Batch Certificate) 포함
- Techno Glass



## Automatic Buret, Pellet Pattern 자동뷰렛, 중간콕크 부착형, Class AS + Batch 보증서

Cat. No.	Model	용량	눈금	EA Price ₩
B13-59-881	22764	10 ml	0.02 ml	445,200
B13-59-885	22767	25 ml	0.05 ml	445,200
B13-59-889	22768	50 ml	0.1 ml	445,200

- Borosilicate Glass 3.3
- PTFE Needle Valve & PTFE Intermediate Stopcock
- 보증서 (Batch Certificate) 포함
- 2 L Buret Bottle 포함
- Brand GmBH



## Automatic Buret, Pellet Pattern 자동뷰렛, Class AS + Batch 보증서

Cat. No.	Model	용량	눈금	EA Price ₩
B13-59-871	22524	10 ml	0.02 ml	355,400
B13-59-875	22528	25 ml	0.05 ml	355,400
B13-59-879	22530	50 ml	0.1 ml	355,400

- Borosilicate Glass 3.3
- PTFE Needle Valve
- 보증서 (Batch Certificate) 포함
- 2 L Buret Bottle 포함
- Brand GmBH

