

Accessory for Automatic Buret 자동 뷰렛 악세서리

Buret Bottle ①

Cat. No.	Model	색상	용량	Joint	EA Price ₩
B13-53-808	436050046	Clear	2 L	29/32	26,300
B13-53-816	436050047	Amber			26,300

- Techno Glass

Rubber Blower ②

Cat. No.	Model	길이	Description	EA Price ₩
B13-53-818	436055002	250 mm	Red Rubber	15,400

Pressure Valve Blower ③

Cat. No.	Model	규격	EA Price ₩
G01-11-037	L1143	Black Rubber	7,600

- Gas Sampling 및 Pumping용
- Pressure 조절용 Valve 부착
- Stainless Steel Ball
- Connector od : Ø6 mm

Double Aspirator Bulb ④

Cat. No.	Model	규격	EA Price ₩
G01-11-911	L0094	Double Aspirator Bulb	9,000

- Automatic Buret 용



Buret Cover ⑤

Cat. No.	Model	재질	Size	EA Price ₩
B13-59-893	16400	PP	id20 x h30 mm	3,000

- Brand GmBH



- Buret Clamp는 253 Page 참조
- Buret Brush는 179 Page 참조
- Buret Stand는 260 Page 참조



Digital Burette / Electric Titrator, dTrite® 디지털 뷰렛 / 적정기, 교반 가능

Cat. No.	Model	용량	최소 분해능	EA Price ₩
B13-238-100	dTrite	99.99 ml	0.01	1,053,300

- 디지털 뷰렛 기능과 교반기 기능을 동시에 수행
- Controller 제어로 반복적인 작업을 자동으로 수행 가능
- Burette를 Controller 제어로 수작업의 오류를 방지할 수 있음
- 우수한 정밀도와 정확성을 제공하는 고품질의 모터 내장
- 분리형 튜빙 및 팁 기본 포함
- DLAB

Cat. No.	B13-238-100
Model	dTrite
Volume Range	0.01mL-99.99mL, Max piston lift is 10mL, resolution is 10μL
Volume Accuracy	R=0.2% CV=0.07%
Velocity	16 stages
Operating Temperature Range	10°C - 30°C
Controller	Burette and Stirrer Program Control

