

Evaporator Adapter for Vial 회전농축기 바이알 어댑터

- 회전 농축기와 바이알을 연결하는 어댑터
- Chemglass

Vial Adapter Set

Cat. No.	Model	Joint	Thread	SET Price ₩
E04-132-040	1318-40	24/40	6ea Set	418,700
E04-132-044	1318-42	29/42		439,200

- 6 종류의 바이알 어댑터들을 모두 포함한 Set

Joint Connecting Adapter Only


Cat. No.	Model	Joint	Thread	EA Price ₩
A01-132-049	1318-10	24/40	20-400	82,800
A01-132-056	1318-12	29/42	(G.P.I.)	93,700

- Glass Part : Outer Joint + Screw G.P.I. 20-400

G.P.I. Thread Connector Only

Cat. No.	Model	Top Thread	Bottom Thread	EA Price ₩
A01-132-608	1318-20	20-400	13-425	42,800
A01-132-616	1318-21		15-425	45,700
A01-132-624	1318-23		20-400	52,500
A01-132-632	1318-24		22-400	53,600
A01-132-640	1318-25		24-400	70,000
A01-132-648	1318-22		18-425	50,700

- Teflon Lined Polypropylene Double Threaded Connector



Rotary Evaporator Adapter for Microwave Vial 회전농축기 바이알 어댑터, 마이크로 웨이브 바이알용

- 구성 : PTFE Joint Adapter, Thread Black Acetal Compression Cap, #017 Viton O-ring
- Microwave Vial은 1245 Page 참조
- Chemglass

PTFE Adapter for Microwave Vial

Cat. No.	Model	Joint	Flange	EA Price ₩
V05-132-590	1318-A-10	24/40	20 mm Crimp	305,100
V05-132-593	1318-A-12	29/42		318,900

Spare Viton O-ring

Cat. No.	Model	Description	PK Price ₩
V05-132-625	305-017	Viton O-ring, #017	12ea 19,000



Inlet Adapter, for Rotary Evaporator 회전농축기용 어댑터

- 회전농축기 사용 전에 건조제나 입자성 물질을 제거하기 위한 별도의 단계가 필요하지 않음
- 스탑콕크 그리스나 유기 잔류물이 접합부로 유출되는 것을 방지하여 접합부 동결 및 손상을 줄일 수 있음
- Chemglass

Cat. No.	Model	Inlet Joint	EA Price ₩
A01-132-201	CG-1327-30	19/38	434,100

Spare Accessory

Cat. No.	Model	Description	EA Price ₩
A01-132-204	CG-1327-31	Glass Inlet Adapter Only	342,700
A01-132-207	CG-1327-32	1/4"-28 Thread Tefzel Union	26,100
A01-132-210	CG-1327-34	1/8" O.D. x 24" Long PTFE Tubing	22,900
A01-132-213	CG-1327-36	20 um Solvent Inlet Filter	58,900
A01-132-216	CG-1327-21	5mm I.D. x 12" Long PTFE Tubing	14,200

