

## Compact Micro Centrifuge, PURISPIN L15 마이크로 원심분리기, 1.5/2.0ml Tube×30개

Cat. No.	Model	구성 Rotor	EA Price
C05-608-031	PURISPIN L15 Set	본체 + A241.5S (Angle Rotor)	3,300,000
C05-608-036	PURISPIN L15	X (Body Only)	2,750,000

- 마이크로 튜브 사용 실험에 적합한 공냉식 원심분리기
- Table Top으로 사용할 수 있는 컴팩트한 디자인
- 0.2ml - 5ml Tube 까지 사용 가능한 Rotor 구성
- Soft Closing Door Lock : 도어 개폐 시 낮은 소음 및 쉽고 부드러운 도어 Locking
- 빠른 가/감속 시간 : 6 Step의 빠른 가/감속 설정
- 작동중에도 모든 파라미터 값 확인 및 변경 가능
- 소음 저감 장치를 통한 저소음 구현
- **CRYTE**

Model	PURISPIN L15
Max. RPM	15,000 RPM
Max. RCF	21,608×g
Max. Capacity	1.5/2.0 ml×30개 / 5 ml×10개
Accel/Decel Phase	0-5 (6 Step)
Noise Level	≤60 dB
Set Time	99 min 59 sec / Continuous
Time Count	From start / From set-up RPM
Program	10 Memory
Display	Blue LCD
RPM Range	400 to 15,000 RPM
Available Rotor Quantity	4 (Angle Rotor)
Display Message on Screen	RPM(RCF), Time, ACC, DEC, Program, Rotor ID
Power Supply	230 V ±10%, 50/60 Hz, 1 Phase / 0.45 KW
Dimension (W×D×H)	291×416×246/550 mm, 19.5 kg (without Rotor)



### Angle Rotor for PURISPIN L15

Cat. No.	C05-608-301	
Model	A480.2S	
Rotor Type	48×0.2 ml	
Max. RPM	11,000 rpm	
Available Tube	0.2 ml	0.2ml 8 Strips×6
Use Adapter	-	
Tube Size	∅6.5 mm	∅6.5 mm
Radius	83.5 mm	83.5 mm
Rcf	11,296×g	11,296×g
Tube Angle	45°	
EA Price ₩	780,000	

Cat. No.	C05-608-307		
Model	A301.5S		
Rotor Type	30×1.5 ml		
Max. RPM	12,000 rpm		
Available Tube	0.2 ml	0.5 ml	1.5/2.0 ml
Use Adapter	F1.5A-0.2	F1.5A-0.5	-
Tube Size	∅6.5×44 mm	∅8×52 mm	∅11.1×61 mm
Radius	80.3 mm	86.7 mm	94.8 mm
Rcf	12,928×g	13,958×g	15,262×g
Tube Angle	45°		
EA Price ₩	780,000		

Cat. No.	C05-608-304		
Model	A241.5S		
Rotor Type	24×1.5 ml		
Max. RPM	15,000 rpm		
Available Tube	0.2 ml	0.5 ml	1.5/2.0 ml
Use Adapter	F1.5A-0.2	F1.5A-0.5	-
Tube Size	∅6.5×44 mm	∅8×51 mm	∅11.1×61 mm
Radius	70.3 mm	76.7 mm	84.8 mm
Rcf	17,684×g	19,294×g	21,331×g
Tube Angle	45°		
EA Price ₩	550,000		

Cat. No.	C05-608-310		
Model	A105S		
Rotor Type	10×5 ml		
Max. RPM	15,000 rpm		
Available Tube	5 ml		
Use Adapter	-		
Tube Size	∅17×67 mm		
Radius	85.9 mm		
Rcf	21,608×g		
Tube Angle	45°		
EA Price ₩	600,000		