

## Refrigerated Micro Centrifuge, PURISPIN L18R 냉장형 마이크로 원심분리기, 50ml Conical × 6개

Cat. No.	Model	구성 Rotor	EA Price
C05-608-021	PURISPIN L18R Set	본체 + A650M (Angle Rotor)	7,400,000
C05-608-026	PURISPIN L18R	X (Body Only)	5,890,000

- 온도에 민감한 샘플도 사용 가능한 냉장형 마이크로 원심분리기
- 0.2ml ~ 50ml까지 사용 가능한 Rotor 구성
- 4.3" Full Color Touch LCD와 터치패드 잠금 기능
- 냉장제어 기술과 로터의 마찰을 최소화하여 정확하고 안정적인 온도 유지 가능
- 비밀번호 설정을 통한 보안 설정
  - 하단의 배수구를 통해 응축수 제거
- AS241.5MX : 병원성 미생물 등 감염 인자가 포함된 에어로졸의 위험을 차단
  - Containment, Certified for dynamic microbiological test (En611010-020) by TUV NORD(독일 글로벌 시험인증기관)
- 각종 Rotor의 자세한 사항은 홈페이지 ([www.labmart.co.kr](http://www.labmart.co.kr)) 참조
- CRYSTE

Model	PURISPIN L18R	
Rotor Type	Angle Rotor	Swing Rotor
Max. RPM	18,000 RPM	13,000 RPM
Max. RCF	30,717×g	16,343×g
Max. Capacity	50 ml Conical×6	1.5/2.0 ml×24
Accel / Decel phase	0 ~ 9 ( 10 Step )	
Set Temp.	-10 ~ 40°C ( 1°C increment )	
Ref. System	Non-CFC	
Noise Level	≤60 dB	
Set Time	99 min 59 sec / Continuous	
Time Count	From start / From set-up RPM	
Program	100 Memory	
Rotor ID	Automatic Recognition	
Display	4.3" LCD Touch Screen	
RPM Range	400 to 18,000 RPM	
Display Message on Screen	RPM, RCF, Time, Door open/close, Temp ACC, DEC, Rotor ID, Program, Comp. On/Off, Key Lock	
Available Rotor Quantity	11 ( Angle / Swing Rotor )	
Power Supply	230 V ± 10%, 50/60 Hz, 1 Phase / 1.5 KW	
Dimension (W×D×H)	382×669×328/684 mm / 52 kg	

### Angle Rotor for PURISPIN L18R

Cat. No.	Model	용량	Max. RPM	EA Price ₩
C05-608-201	A480.2M	0.2 ml×48개	11,000 rpm	820,000
C05-608-204	A241.5M	1.5/2.0 ml×24개	18,000 rpm	650,000
C05-608-207	AS241.5MX	1.5/2.0 ml×24개 (Aerosol-tight)		1,340,000
C05-608-210	A301.5M	1.5/2.0 ml×30개	14,000 rpm	820,000
C05-608-222	A125M	5 ml×12개		680,000
C05-608-219	A181.5MH	1.5/2.0 ml×18개	15,000 rpm	910,000
C05-608-213	A481.5M	1.5/2.0 ml×48개	13,000 rpm	960,000
C05-608-225	A650M	50 ml×6개	8,000 rpm	1,510,000

### Swing Rotor for PURISPIN L18R

Cat. No.	Model	용량	Wing	Max. RPM	EA Price ₩
C05-608-228	S61.5M	1.5/2.0 ml×24개	6개	13,000 rpm	1,490,000
C05-608-231	S2M-MTP	Microplate용	2개	5,000 rpm	2,520,000

### Swing Rotor Bucket for PURISPIN L18R

Cat. No.	Model	용량	적용 Rotor	PK Price ₩
C05-608-558	B1.5	1.5/2.0 ml×4개	S61.5M	6ea 180,000

### Tube Adapter

Cat. No.	Model	용량	적용 Rotor	PK Price ₩
C05-608-501	F1.5A-0.2	0.2 ml	A241.5M, AS241.5MX, A301.5M, A181.5MH, A481.5M, S61.5M	2ea 6,000
C05-608-504	F1.5A-0.5	0.5 ml	A241.5M, AS241.5MX, A301.5M, A181.5MH, A481.5M, S61.5M	2ea 6,000
C05-608-510	F50CA-15C	15 ml Conical	A650M	2ea 60,000
C05-608-605	S50A-50C	50 ml Conical		2ea 24,000

