

High Speed Centrifuge, VELOSPIN L22R 대용량 고속 원심분리기, Max. 1L×6개

Cat. No.	Model	용량	Max. RPM	EA Price
C05-608-011	VELOSPIN L22R	1L × 6개	22,000RPM	17,000,000

- 강력한 파워와 대용량 시료 처리가 가능한 냉장형 고속 원심분리기
- 0.2 ml Tube에서 1L Bottle까지 다양한 용기 사용 가능
- 4.3" Full Color Touch LCD와 터치패드 잠금 기능
- 냉장제어 기술과 로터의 마찰을 최소화하여 정확하고 안정적인 온도 유지 가능
- Rotor의 미끄러짐 방지를 위한 안전 허브 시스템으로 높은 안전성 유지
- 저소음, 냉장 및 신속 냉각 • 비밀번호 설정을 통한 보안 설정
- 하단의 배수구를 통해 응축수 제거
- 각종 Rotor의 자세한 사항은 홈페이지 (www.labmart.co.kr) 참조
- **CRYTE**

Model	VELOSPIN L22R	
Rotor Type	Angle Rotor	Swing Rotor
Max. RPM	22,000 RPM	4,000 RPM
Max. RCF	54,057×g	3,265×g
Max. Capacity	1,000 ml×6	250 ml×4
Accel / Decel phase	0 ~ 9 (10 Step)	
Set Temp.	-10 ~ 40°C (1°C increment)	
Ref. System	Non-CFC	
Noise Level	≤60 dB	
Set Time	99 hour 59 min 59 sec / Continuous	
Time Count	From start / From set-up RPM	
Program	100 Memory	
Rotor ID	Automatic Recognition	
Display	4.3" LCD Touch Screen	
RPM Range	400 to 22,000 RPM	
Display Message on Screen	RPM, RCF, Time, Door open/close, Temp ACC, DEC, Rotor ID, Comp. On/Off, Key Lock	
Available Rotor Quantity	14 (Angle / Swing Rotor)	
Power Supply	230 V ± 10%, 50/60 Hz, 1 Phase / 5.8 KW	
Dimension (W×D×H)	w667×d772×h890/1530 mm / 232 kg	

Angle Rotor for VELOSPIN L22R

Cat. No.	Model	용량	Max. RPM	EA Price ₩
C05-608-101	A480.2H	0.2 ml×48개	12,000 rpm	1,210,000
C05-608-104	A301.5H	1.5/2.0 ml×30개	22,000 rpm	1,210,000
C05-608-107	A481.5H	1.5/2.0 ml×48개	13,000 rpm	1,330,000
C05-608-110	A1215H	15 ml×12개	20,000 rpm	2,420,000
C05-608-113	A1215HC	15 ml Conical×12개	20,000 rpm	2,420,000
C05-608-116	A650H	50 ml×6개	22,000 rpm	2,450,000
C05-608-119	A650HC	50 ml Conical×6개	17,000 rpm	2,470,000
C05-608-122	A850H	50 ml×8개	20,000 rpm	2,860,000
C05-608-125	A685H	85 ml×6개	20,000 rpm	2,640,000
C05-608-128	A6250H	250 ml×6개	12,000 rpm	4,400,000
C05-608-131	A6500H	500 ml×6개	10,000 rpm	5,900,000
C05-608-134	A41000H	1000 ml×4개	8,000 rpm	8,740,000
C05-608-137	A61000H	1000 ml×6개	7,000 rpm	10,890,000

Swing Rotor & Bucket for VELOSPIN L22R

Cat. No.	Model	최대 용량	Wing	Max. RPM	EA Price ₩
C05-608-140	S4250HE	250 ml×4개	4개	4,000 rpm	980,000

Tuberack for "EB250" & "EBS250" Bucket

Cat. No.	Model	용량	Adapter Bore	PK Price ₩
C05-608-531	S250ET-3-12	3 ml×12개	Ø13.2 mm	2ea 100,000
C05-608-534	S250ET-5-12	5 ml×12개		2ea 100,000
C05-608-537	S250ET-10-12	10 ml×12개	Ø16 mm	2ea 100,000
C05-608-540	S250ET-15-12	15 ml×12개	Ø17.2 mm	2ea 100,000
C05-608-543	S250ET-15C-9	15 ml Conical×9개	Ø17 mm	2ea 100,000
C05-608-546	S250ET-50-4	50 ml×4개		2ea 100,000
C05-608-549	S250ET-50C-3	50 ml Conical×3개	Ø29.5 mm	2ea 100,000
C05-608-552	S250ET-85-2	85 ml×2개	Ø38.5 mm	2ea 100,000
C05-608-555	S250EA-250	250 ml×1개	Ø62.5 mm	2ea 100,000

• Tube Adapter는 224 ~ 225 Page 참조



Bucket for "S4250HE" Swing Out Rotor

Cat. No.	Model	Type	Max. RPM	PK Price ₩
C05-608-522	EB250	Rectangular	4,000 rpm	4ea 960,000
C05-608-525	EBS250	Biocontainment		4ea 1,520,000
C05-608-528	BW250	Microplate		4ea 800,000