

## Refrigerated Multi-Purpose Centrifuge, VARISPIN L15R 중형 다목적 냉장형 원심분리기, 750 ml×6개

Cat. No.	Model	구성품	EA Price
C05-608-041	VARISPIN L15R Set	Rotor, Bucket (4개), Tuberack (4개)	9,900,000
C05-608-046	VARISPIN L15R	Body Only	7,150,000

- 강력한 파워와 폭넓은 사용이 가능한 다목적 냉장 원심분리기
- 0.2 ml Tube에서 750ml Bottle까지 다양한 용기 사용 가능
- 4.3" Full Color Touch LCD와 터치패드 잠금 기능
- "EBS250" & "BS750" Bucket :  
병원성 미생물 등 감염 인자가 포함된 에어로졸의 위험을 차단  
- Containment, Certified for dynamic microbiological test (En61010-020)  
by TUV NORD(독일 글로벌 시험인증기관)
- Set 구성품 : Main Body, Swing Rotor (S4750T) 1개, Bucket (B750) 4개,  
Tuberack 4개로 구성
- 각종 Rotor의 자세한 사항은 홈페이지 (www.labmart.co.kr) 참조
- CRYSTE

Model	VARISPIN L15R	
Rotor Type	Angle Rotor	Swing Rotor
Max. RPM	15,000 RPM	5,000 RPM
Max. RCF	23,847×g	4,947×g
Max. Capacity	250 ml×6	750 ml×4
Accel / Decel phase	0 ~ 9 ( 10 Step )	
Set Temp.	-10 ~ 40°C ( 1°C increment )	
Ref. System	Non-CFC	
Noise Level	≤60 dB	
Set Time	99 hour 59 min 59 sec / Continuous	
Time Count	From start / From set-up RPM	
Program	100 Memory	
Rotor ID	Automatic Recognition	
Display	4.3" LCD Touch Screen	
RPM Range	400 to 15,000 RPM	
Display Message on Screen	RPM, RCF, Time, Door open/close, Temp ACC, DEC, Rotor ID, Program, Comp. On/Off, Key Lock	
Available Rotor Quantity	16 ( Angle / Swing Rotor )	
Power Supply	230 V ± 10%, 50/60 Hz, 1 Phase / 2.5 KW	
Dimension (W×D×H)	780×687×375/915 mm / 98 kg	

### Angle Rotor for VARISPIN L15R

Cat. No.	Model	용량	Max. RPM	EA Price ₩
C05-608-401	A480.2T	0.2 ml×48개	12,000 rpm	1,160,000
C05-608-404	A301.5T	1.5/2.0 ml×30개	15,000 rpm	1,160,000
C05-608-407	A481.5T	1.5/2.0 ml×48개	13,000 rpm	1,270,000
C05-608-410	A1215T	15 ml×12개	12,000 rpm	2,310,000
C05-608-413	A1215TC	15 ml Conical×12개	13,000 rpm	2,310,000
C05-608-416	A650T	50 ml×6개	12,000 rpm	2,180,000
C05-608-419	A650TC	50 ml Conical×6개	13,000 rpm	2,180,000
C05-608-422	A850T	50 ml×8개	12,000 rpm	2,310,000
C05-608-425	A685T	85 ml×6개	12,000 rpm	2,310,000
C05-608-450	A4250T	250 ml×4개	8,000 rpm	2,600,000
C05-608-428	AS301.5TX	1.5/2.0 ml×30개 (Aerosol-tight)	15,000 rpm	2,100,000
C05-608-431	AS685TX	85 ml×6개 (Aerosol-tight)	12,000 rpm	4,000,000

### Swing Rotor only for VARISPIN L15R

Cat. No.	Model	최대 용량	Wing	Bucket(별매)	EA Price ₩
C05-608-434	S4250TE	250 ml×4개	4개	EB250, EBS250, BW250	970,000
C05-608-437	S4250T	250 ml×4개	4개	B250	1,110,000
C05-608-440	S4750T	750 ml×4개	4개	B750, BS750, B500, BW750	1,710,000
C05-608-443	S4100T	100 ml×4개	4개	B100	1,080,000

- Swing Rotor 만의 가격이며 Bucket 및 Tube Adapter는 별도 구매해야 함



AS301.5TX



AS685TX



S4750T



S4100T



B100

### Bucket for Swing Rotor

Cat. No.	Model	적용 Rotor	Type	PK Price ₩
C05-608-522	EB250	S4250TE	Rectangular	4ea 960,000
C05-608-525	EBS250		Biocontainment	4ea 1,520,000
C05-608-528	BW250	S4250T	Microplate	4ea 800,000
C05-608-561	B250		Round	4ea 680,000
C05-608-564	B750	S4750T	Round	4ea 800,000
C05-608-567	BS750		Biocontainment	4ea 2,000,000
C05-608-570	B500		Round	4ea 1,200,000
C05-608-573	BW750		Microplate	4ea 1,280,000
C05-608-576	B100	S4100T	Round	4ea 520,000

### Swing Rotor 별 사용 Tube (용기)

Swing Rotor	Bucket	사용 Tube
S4250TE	EB250, EBS250	3 ml, 5 ml, 10 ml, 15 ml, 15 ml Conical, 50 ml, 50 ml Conical, 85 ml, 250 ml
	BW250	Microplate (MTP, DWP)
S4250T	B250	5 ml, 10 ml, 15 ml, 15 ml Conical, 50 ml, 50 ml Conical, 250 ml
S4750T	B750, BS750	3 ml, 5 ml, 10 ml, 15 ml, 15 ml Conical, 50 ml, 50 ml Conical, 250 ml, 750 ml
	B500	500 ml
	BW750	Microplate (MTP, DWP)
S4100T	B100	1.5 ml, 5 ml, 10 ml, 15 ml, 15 ml Conical, 50 ml, 50 ml Conical, 85 ml, 100 ml

- Tube Adapter는 224 ~ 225 Page 참조