

Refrigerated Multi-Purpose Centrifuge, VARISPIN L12R 중형 다목적 냉장형 원심분리기, 250 ml×4개

Cat. No.	Model	구성품	EA Price
C05-608-061	VARISPIN L12R Set	Rotor, Bucket (4개), Tuberack (4개)	8,700,000
C05-608-066	VARISPIN L12R	Body Only	6,570,000

- 많은 양의 시료 처리와 폭넓은 사용이 가능한 다목적 냉장 원심분리기
- 냉장제어 기술과 로터의 마찰을 최소화하여 정확하고 안정적인 온도 유지 가능
- 0.2 ml Tube에서 250ml Bottle까지 다양한 용기 사용 가능
- 4.3" Full Color Touch LCD와 터치패드 잠금 기능
- "EBS250" Bucket : 병원성 미생물 등 감염 인자가 포함된 에어로졸의 위험을 차단
- Containment, Certified for dynamic microbiological test (En61010-020) by TUV NORD(독일 글로벌 시험인증기관)
- Set 구성품 : Main Body, Swing Rotor (S4250TE) 1개, Bucket (EB250) 4개, Tuberack (S250ET-50C-3×2개, S250ET-15C-9×2개) 4개로 구성
- 각종 Rotor의 자세한 사항은 홈페이지 (www.labmart.co.kr) 참조
- **CRYTE**

Model	VARISPIN L12R	
Rotor Type	Angle Rotor	Swing Rotor
Max. RPM	12,000 RPM	5,000 RPM
Max. RCF	16,647 × g	4,947 × g
Max. Capacity	85 ml×6	250 ml×4
Accel / Decel phase	0 ~ 9 (10 Step)	
Set Temp.	-10 ~ 40°C (1°C increment)	
Ref. System	Non-CFC	
Noise Level	≤60 dB	
Set Time	99 hour 59 min 59 sec / Continuous	
Time Count	From start / From set-up RPM	
Program	100 Memory	
Rotor ID	Automatic Recognition	
Display	4.3" LCD Touch Screen	
RPM Range	400 to 12,000 RPM	
Display Message on Screen	RPM, RCF, Time, Door open/close, Temp ACC, DEC, Rotor ID, Program, Comp. On/Off, Key Lock	
Available Rotor Quantity	14 (Angle / Swing Rotor)	
Power Supply	230 V ± 10%, 50/60 Hz, 1 Phase / 2 KW	
Dimension (W×D×H)	700×642×375/895 mm / 91 kg	

Angle Rotor for VARISPIN L12R

Cat. No.	Model	용량	Max. RPM	EA Price ₩
C05-608-401	A480.2T	0.2 ml×48개	12,000 rpm	1,160,000
C05-608-404	A301.5T	1.5/2.0 ml×30개		1,160,000
C05-608-407	A481.5T	1.5/2.0 ml×48개		1,270,000
C05-608-410	A1215T	15 ml×12개		2,310,000
C05-608-413	A1215TC	15 ml Conical×12개		2,310,000
C05-608-416	A650T	50 ml×6개		2,180,000
C05-608-419	A650TC	50 ml Conical×6개		2,180,000
C05-608-422	A850T	50 ml×8개		2,310,000
C05-608-425	A685T	85 ml×6개		2,310,000
C05-608-428	AS301.5TX	1.5/2.0 ml×30개 (Aerosol-tight)		2,100,000
C05-608-431	AS685TX	85 ml×6개 (Aerosol-tight)		4,000,000

Swing Rotor only for VARISPIN L12R

Cat. No.	Model	최대 용량	Wing	Bucket(별매)	EA Price ₩
C05-608-434	S4250TE	250 ml×4개	4개	EB250, EBS250, BW250	970,000
C05-608-437	S4250T	250 ml×4개	4개	B250	1,110,000
C05-608-443	S4100T	100 ml×4개	4개	B100	1,080,000

• Swing Rotor 만의 가격이며 Bucket 및 Tube Adapter는 별도 구매해야 함



Bucket for Swing Rotor

Cat. No.	Model	적용 Rotor	Type	PK Price ₩
C05-608-522	EB250	S4250TE	Rectangular	4ea 960,000
C05-608-525	EBS250		Biocontainment	4ea 1,520,000
C05-608-528	BW250		Microplate	4ea 800,000
C05-608-561	B250	S4250T	Round	4ea 680,000
C05-608-576	B100	S4100T	Round	4ea 520,000

Swing Rotor 별 사용 Tube (용기)

Swing Rotor	Bucket	사용 Tube
S4250TE	EB250, EBS250	3 ml, 5 ml, 10 ml, 15 ml, 15 ml Conical, 50 ml, 50 ml Conical, 85 ml, 250 ml
	BW250	Microplate (MTP, DWP)
S4250T	B250	5 ml, 10 ml, 15 ml, 15 ml Conical, 50 ml, 50 ml Conical, 250 ml
S4100T	B100	1.5 ml, 5 ml, 10 ml, 15 ml, 15 ml Conical, 50 ml, 50 ml Conical, 85 ml, 100 ml

• Tube Adapter는 224 ~ 225 Page 참조