

Refrigerated Multi-Purpose Centrifuge, VARISPIN L6R 중형 다목적 냉장형 원심분리기, 100 ml×4개

Cat. No.	Model	구성품	EA Price
C05-608-761	VARISPIN L6R Set	Rotor, Bucket (4개), Tuberack (4개)	7,000,000
C05-608-764	VARISPIN L6R	Body Only	5,640,000

- 중급 용량의 시료 처리와 폭넓은 사용이 가능한 다목적 냉장 원심분리기
- 냉장제어 기술과 로터의 마찰을 최소화하여 정확하고 안정적인 온도 유지 가능
- 1.5 ml Tube에서 100ml Bottle까지 다양한 용기 사용 가능
- 4.3" Full Color Touch LCD
- Set 구성품 : Main Body, Swing Rotor (S4100G) 1개, Bucket (B100) 4개, Tuberack 4개로 구성
- 각종 Rotor의 자세한 사항은 홈페이지 (www.labmart.co.kr) 참조
- **CRYSTE**



Model	VARISPIN L6R	
Rotor Type	Angle Rotor	Swing Rotor
Max. RPM	6,000 RPM	4,500 RPM
Max. RCF	5,357×g	3,781×g
Max. Capacity	50 ml×8	100 ml×4
Accel / Decel phase	0 ~ 9 (10 Step)	
Set Temp.	-10 ~ 40°C (1°C increment)	
Ref. System	Non-CFC	
Noise Level	≤60 dB	
Set Time	99 hour 59 min 59 sec / Continuous	
Time Count	From start / From set-up RPM	
Program	100 Memory	
Rotor ID	Automatic Recognition	
Display	4.3" LCD Touch Screen	
RPM Range	400 to 6,000 RPM	
Display Message on Screen	RPM, RCF, Time, Door open/close, Temp ACC, DEC, Rotor ID, Program, Comp. On/Off, Key Lock	
Available Rotor Quantity	4 (Angle / Swing Rotor)	
Power Supply	230 V ±10%, 50/60 Hz, 1 Phase / 1.8 KW	
Dimension (W×D×H)	636×611×334/774 mm / 67 kg	

Angle Rotor for VARISPIN L6R

Cat. No.	Model	용량	Max. RPM	EA Price ₩
C05-608-610	AL2415G	15 ml×24개	4,000 rpm	550,000
C05-608-640	A850GC	50 ml Conical×8개	6,000 rpm	1,730,000

Swing Rotor only for VARISPIN L6R

Cat. No.	Model	최대 용량	Wing	Bucket(별매)	EA Price ₩
C05-608-620	S4100G	100 ml×4개	4개	B100	640,000
C05-608-628	S2G-MTP	Microplate×2개	2개	-	460,000

- Swing Rotor 만의 가격이며 Bucket 및 Tube Adapter는 별도 구매해야 함

Bucket for Swing Rotor

Cat. No.	Model	적용 Rotor	Type	PK Price ₩
C05-608-576	B100	S4100G	Round	4ea 520,000

Swing Rotor 별 사용 Tube (용기)

Swing Rotor	Bucket	사용 Tube
S4100G	B100	1.5 ml, 5 ml, 10 ml, 15 ml, 15 ml Conical, 50 ml, 50 ml Conical, 85 ml, 100 ml

- Tube Adapter는 224 ~ 225 Page 참조

