

Multi Centrifuge, VARISPIN L6 다목적 공냉식 원심분리기, 250 ml×4개

Cat. No.	Model	구성품	EA Price
C05-608-096	VARISPIN L6 Set	Rotor, Bucket (4개), Tuberack (4개)	4,700,000
C05-608-098	VARISPIN L6	Body Only	3,340,000

- 컴팩트한 디자인에 높은 시료 처리 능력을 가짐
- Set 구성품 : Main Body, Swing Rotor (S4100G) 1개, Bucket (B100) 4개, Tuberack 4개로 구성
- 각종 Rotor의 자세한 사항은 홈페이지 (www.labmart.co.kr) 참조
- CRYSTE

Model	VARISPIN L6	
Rotor Type	Angle Rotor	Swing Rotor
Max. RPM	6,000 RPM	5,000 RPM
Max. RCF	5,357×g	4,724×g
Max. Capacity	50 ml×8	250 ml×4
Accel / Decel phase	0 ~ 5 (6 Step)	
Noise Level	≤60 dB	
Set Time	99 min 59 sec / Continuous	
Time Count	From start / From set-up RPM	
Program	10 Memory	
Display	Blue LCD	
RPM Range	400 to 6,000 RPM	
Display Message on Screen	RPM, RCF, Time, ACC, DEC, Program, Rotor ID	
Available Rotor Quantity	7 (Angle / Swing Rotor)	
Power Supply	230 V ±10%, 50/60 Hz, 1 Phase / 1.1 KW	
Dimension (W×D×H)	427×560×320/731 mm / 32 kg	

Angle Rotor for VARISPIN L6

Cat. No.	Model	용량	Max. RPM	EA Price ₩
C05-608-615	AL1615G	15 ml×16개	4,000 rpm	480,000
C05-608-610	AL2415G	15 ml×24개		550,000
C05-608-640	A850GC	50 ml×8개	6,000 rpm	1,730,000

- AL1615G & AL2415G : 15 ml용 Stainless Steel Sleeve (Ø18 mm) 포함
- 3 ml, 5 ml, 10 ml, 15 ml, 15 ml Conical Tube 용

Tube Adapter for Angle Rotor

Cat. No.	Model	전용 Rotor	용량	Bore	PK Price ₩
C05-608-634	F15A-3	AL1615G, AL2415G	3 ml 용	13.5 mm	2ea 12,000
C05-608-636	F15A-5		5 ml 용		2ea 12,000
C05-608-510	F50CA-15C	A850GC	15 ml 용	17 mm	2ea 60,000

Swing Rotor for VARISPIN L6

Cat. No.	Model	용량	Wing	Max. RPM	EA Price ₩
C05-608-620	S4100G	100 ml×4개	4개	4,000 rpm	640,000
C05-608-625	S650G	50 ml×6개	6개		640,000
C05-608-644	S4250GE	250 ml×4개	4개	3,500 rpm	730,000
C05-608-628	S2G-MTP	Microplate×2개	2개	4,000 rpm	460,000

- S4100G & S650G Rotor 및 Tube Adapter 223 ~ 225 Page 참조
- S4250GE : 3 ml, 5 ml, 10 ml, 15 ml, 15 ml Conical,
50 ml, 50ml Conical, 85 ml, 250 ml용

Tuberack for Swing Rotor "S4250GE"

Cat. No.	Model	용량	Adapter Bore	PK Price ₩
C05-608-531	S250ET-3-12	3 ml×12개	Ø13.2 mm	2ea 100,000
C05-608-534	S250ET-5-12	5 ml×12개		2ea 100,000
C05-608-537	S250ET-10-12	10 ml×12개	Ø16 mm	2ea 100,000
C05-608-540	S250ET-15-12	15 ml×12개	Ø17.2 mm	2ea 100,000
C05-608-545	S250EG-15C-7	15 ml Conical×7개	Ø17 mm	2ea 360,000
C05-608-546	S250ET-50-4	50 ml×4개	Ø29.5 mm	2ea 100,000
C05-608-548	S250EG-50C-3	50 ml Conical×3개		2ea 340,000
C05-608-552	S250ET-85-2	85 ml×2개	Ø38.5 mm	2ea 100,000
C05-608-555	S250EA-250	250 ml×1개	Ø62.5 mm	2ea 100,000



Bucket for Swing Rotor

Cat. No.	Model	적용 Rotor	Type	PK Price ₩
C05-608-576	B100	S4100G	Round	4ea 520,000
C05-608-646	B50	S650G	Round	6ea 420,000
C05-608-522	EB250	S4250GE	Rectangular	4ea 960,000