


CRYSTE Tube Adapter for Rotor 에서 계속...
"B250" Round Bucket 용 Tuberack

Cat. No.	Model	용량	Bore	Max. Height	PK Price ₩
C05-608-728	S250T-5-8	5 ml×8개	Ø13.5mm	145 mm	2ea 140,000
C05-608-731	S250T-10-7	10 ml×7개	Ø16mm		2ea 140,000
C05-608-734	S250T-15-5	15 ml×5개	Ø17.5mm		2ea 140,000
C05-608-737	S250T-15C-4	15 ml Conical ×4개	Ø17mm	148 mm	2ea 140,000
C05-608-740	S250T-50-2	50 ml×2개	Ø29.5mm	145 mm	2ea 140,000
C05-608-743	S250A-50C	50 ml Conical ×1개	Ø30mm	150 mm	2ea 140,000


"B50" Round Bucket 용 Tuberack

Cat. No.	Model	용량	Bore	Max. Height	PK Price ₩
C05-608-649	S50A-15C	15 ml Conical ×1개	Ø17mm	126 mm	2ea 24,000
C05-608-652	S50A-15	15 ml×1개	Ø17.5mm	123 mm	2ea 24,000
C05-608-605	S50A-50C	50 ml Conical ×1개	Ø30mm	126 mm	2ea 24,000


**Clinical Master Centrifuge, Pro-Cyt
임상 검사용 원심분리기, 실험실용**

Cat. No.	Model	Type	EA Price ₩
C05-79-195	PrO-Cyt LCD	Cytology Centrifuge	6,184,600

- LCD Display
- 1 rpm 단위 속도 조절
- 55dB 이하의 소음
- Emergency Lid Release
- Rotor 발매
- 10 Memory Program
- 0 ~ 99 min 타이머 (30초 단위 조절)
- Multi Point Lid Locking
- CE / ISO9001 / IVD Certified
- **Centurion Scientific Ltd.**

Cat. No.	C05-79-195
Model	PrO-Cyt LCD
Size	w325×d470×h275 mm
Speed	200 ~ 2,000 rpm
Max. Rcf	550 G
Power	175 watts


Optional Rotor, Autoclavable


Cat. No.	C05-79-669
Model	4420
Rotor Type	4 x 0.2 ~ 6ml
Holders	Single or Double
Speed	200 ~ 2,000 rpm
Max. Rcf	550 G
Radius Max	12
Tube Angle	0°
EA Price ₩	3,019,000



Cat. No.	C05-79-671
Model	4430
Rotor Type	8 x 0.2 ~ 6ml
Holders	Single or Double
Speed	200 ~ 2,000 rpm
Max. Rcf	550 G
Radius Max	12
Tube Angle	0°
EA Price ₩	4,006,900



Cat. No.	C05-79-673
Model	4460
Rotor Type	12 x 0.2 ~ 6ml
Holders	Single or Double
Speed	200 ~ 2,000 rpm
Max. Rcf	550 G
Radius Max	12
Tube Angle	0°
EA Price ₩	5,344,200

Optional Accessory

Cat. No.	Model	Type	PK Price ₩
C05-79-980	4444	Single Sample Holder with Card (up to 1 ml)	24ea 185,200
C05-79-982	4446	Double Sample Holder with Card (up to 1 ml)	24ea 220,200
C05-79-984	4600	Single Sample Holder with Card (up to 6 ml)	40ea 595,500
C05-79-986	4700	Stainless Steel Clip	1ea 149,400

