

Micro Centrifuge, Pro-Analytical 마이크로 원심분리기, Max. 120ml

Cat. No.	Model	Type	EA Price ₩
C05-79-101	C1006	Micro Dual Centrifuge	2,000,700
C05-79-116	C1015	Micro Prime Centrifuge	2,325,900

- LED Display
- 10 RPM 단위 속도조절
- 60 db 이하의 소음
- Emergency Lid-Release
- Rotor 별매
- 10 Memory Program
- 0 ~ 99 min 타이머 (30초 단위 조절가능)
- Multi-Point Lid-Locking
- CE / ISO9001 / IVD Certified
- Centurion Scientific Ltd.

Cat. No.	C05-79-101	C05-79-116
Model	C1006	C1015
Size	w235×d350×h235 mm / 12.8 kg	
Speed	500 ~ 8,000 rpm	500 ~ 15,000 rpm
Max. Rcf	6,100 G	22,000 G
Power	140 watts	



Reducer는 230 Page 참조



Optional Rotor, Autoclavable at 121°C



Cat. No.	C05-79-643
Model	BRK 5407
Rotor Type	4×15 ml & 4×10 ml
Speed	500 ~ 8000 rpm
Tube Size Max.	17×120 / 13×110 mm
Max. Rcf	6,100 G
Radius Max	8.5
Tube Angle	30°
Autoclavable	121°C
사용 원심분리기	C1006 / C1015
EA Price ₩	693,100



Cat. No.	C05-79-645
Model	BRK 5408
Rotor Type	8×15 ml
Speed	500 ~ 6000 rpm
Tube Size Max.	17×120 mm
Max. Rcf	3,500 G
Radius Max	8.5
Tube Angle	30°
Autoclavable	121°C
사용 원심분리기	C1006 / C1015
EA Price ₩	762,400



Cat. No.	C05-79-545
Model	BRK 5401
Rotor Type	24 Capillary & 12×2 ml
Speed	500 ~ 12,000 rpm
Tube Size Max.	2×75 / 11×40 mm
Max. Rcf	13,500 G
Radius Max	8.5
Tube Angle	0 & 60°
Autoclavable	121°C
사용 원심분리기	C1015
EA Price ₩	930,800



Cat. No.	C05-79-549
Model	BRK 5424
Rotor Type	24×2 ml
Speed	500 ~ 15,000 rpm
Tube Size Max.	11×50 mm
Max. Rcf	22,000 G
Radius Max	8.5
Tube Angle	45°
Autoclavable	121°C
사용 원심분리기	C1015
EA Price ₩	877,200



Cat. No.	C05-79-653
Model	BRK 5436
Rotor Type	36×0.5 ml
Speed	500 ~ 15,000 rpm
Tube Size Max.	8×30 mm
Max. Rcf	22,000 G
Radius Max	8.5
Tube Angle	45°
Autoclavable	121°C
사용 원심분리기	C1015
EA Price ₩	971,200



Cat. No.	C05-79-553
Model	BRK 5448
Rotor Type	48×0.2 ml
Speed	500 ~ 15,000 rpm
Tube Size Max.	6×40 mm
Max. Rcf	22,000 G
Radius Max	8.5
Tube Angle	45°
Autoclavable	121°C
사용 원심분리기	C1015
EA Price ₩	930,800



Cat. No.	C05-79-655
Model	BRK 5494
Rotor Type	4×PCR Strips
Speed	500 ~ 15,000 rpm
Tube Size Max.	6×40 mm
Max. Rcf	22,000 G
Radius Max	8.5
Tube Angle	45°
Autoclavable	121°C
사용 원심분리기	C1015
EA Price ₩	971,200



Cat. No.	C05-79-675
Model	BRK 5508S
Rotor Type	8×5 ml Swing Out
Speed	500 ~ 4,000 rpm
Tube Size Max.	13×17 mm
Max. Rcf	1,500 G
Radius Max	8.5
Tube Angle	0°
Autoclavable	121°C
사용 원심분리기	C1006 / C1015
EA Price ₩	693,100