

## Small Centrifuge, Pro-Analytical 소형 원심분리기, Max. 400ml

| Cat. No.   | Model  | Type                     | EA Price ₩ |
|------------|--------|--------------------------|------------|
| C05-79-121 | C2004  | Small Centrifuge         | 2,402,000  |
| C05-79-140 | CR2000 | Multi Purpose Centrifuge | 2,953,500  |

- LED Display
- 10 RPM 단위 속도조절
- 60 db 이하의 소음
- Emergency Lid-Release
- Rotor 별매
- 10 Memory Program
- 0 ~ 99 min 타이머 (30초 단위 조절가능)
- Multi-Point Lid-Locking
- CE / ISO9001 / IVD Certified
- Centurion Scientific Ltd.

| Cat. No. | C05-79-121                 | C05-79-140       |
|----------|----------------------------|------------------|
| Model    | C2004                      | CR2000           |
| Size     | w325×d470×h275 mm/ 21.8 kg |                  |
| Speed    | 500 ~ 4,000 rpm            | 500 ~ 15,000 rpm |
| Max. Rcf | 2,200 G                    | 22,200 G         |
| Power    | 160 watts                  | 175 watts        |



CE

Reducer는 230 Page 참조



### Optional Rotor, Autoclavable at 121°C



| Cat. No.       | C05-79-639      |
|----------------|-----------------|
| Model          | BRK 5324        |
| Rotor Type     | 24×15 ml        |
| Speed          | 500 ~ 4,000 rpm |
| Tube Size Max. | 17×120 mm       |
| Max. Rcf       | 2,200 G         |
| Radius Max     | 12.3            |
| Tube Angle     | 30°             |
| Autoclavable   | 121°C           |
| 사용 원심분리기       | C2004 / CR2000  |
| EA Price ₩     | 1,274,200       |



| Cat. No.       | C05-79-637      |
|----------------|-----------------|
| Model          | BRK 5308        |
| Rotor Type     | 8×50 ml         |
| Speed          | 500 ~ 4,000 rpm |
| Tube Size Max. | 30×120 mm       |
| Max. Rcf       | 2,200 G         |
| Radius Max     | 12.3            |
| Tube Angle     | 30°             |
| Autoclavable   | 121°C           |
| 사용 원심분리기       | C2004 / CR2000  |
| EA Price ₩     | 1,228,100       |



| Cat. No.        | C05-79-663      |
|-----------------|-----------------|
| Model           | BRK 1010        |
| Rotor Type      | 4×100 ml        |
| Speed           | 500 ~ 4,000 rpm |
| B5010 Bucket 포함 |                 |
| Max. Rcf        | 2,200 G         |
| Radius Max      | 12.3            |
| Tube Angle      | 0°              |
| Autoclavable    | 121°C           |
| 사용 원심분리기        | C2004           |
| EA Price ₩      | 1,861,700       |



| Cat. No.       | C05-79-647     |
|----------------|----------------|
| Model          | BRK 5412       |
| Rotor Type     | 12×15 ml       |
| Speed          | 500 ~ 6000 rpm |
| Tube Size Max. | 17×120 mm      |
| Max. Rcf       | 3,800 G        |
| Radius Max     | 9.5            |
| Tube Angle     | 30°            |
| Autoclavable   | 121°C          |
| 사용 원심분리기       | C2006          |
| EA Price ₩     | 999,300        |



| Cat. No.       | C05-79-641     |
|----------------|----------------|
| Model          | BRK 5406       |
| Rotor Type     | 6×50 ml        |
| Speed          | 500 ~ 6000 rpm |
| Tube Size Max. | 30×120 mm      |
| Max. Rcf       | 3,800 G        |
| Radius Max     | 9.5            |
| Tube Angle     | 30°            |
| Autoclavable   | 121°C          |
| 사용 원심분리기       | C2006          |
| EA Price ₩     | 890,500        |



| Cat. No.       | C05-79-661     |
|----------------|----------------|
| Model          | BRK 5508M      |
| Rotor Type     | 8×10 ml        |
| Speed          | 500 ~ 4000 rpm |
| Tube Size Max. | 15×115 mm      |
| Max. Rcf       | 2,200 G        |
| Radius Max     | 12.3           |
| Tube Angle     | 0°             |
| Autoclavable   | 121°C          |
| 사용 원심분리기       | CR2000         |
| EA Price ₩     | 921,100        |



| Cat. No.       | C05-79-545             |
|----------------|------------------------|
| Model          | BRK 5401               |
| Rotor Type     | 24 Capillary & 12×2 ml |
| Speed          | 500 ~ 12,000 rpm       |
| Tube Size Max. | 2×75 / 11×40 mm        |
| Max. Rcf       | 13,500 G               |
| Radius Max     | 8.5                    |
| Tube Angle     | 0 & 60°                |
| Autoclavable   | 121°C                  |
| 사용 원심분리기       | CR2000                 |
| EA Price ₩     | 930,800                |



| Cat. No.       | C05-79-549       |
|----------------|------------------|
| Model          | BRK 5424         |
| Rotor Type     | 24×2 ml          |
| Speed          | 500 ~ 15,000 rpm |
| Tube Size Max. | 11×50 mm         |
| Max. Rcf       | 22,000 G         |
| Radius Max     | 8.5              |
| Tube Angle     | 45°              |
| Autoclavable   | 121°C            |
| 사용 원심분리기       | CR2000           |
| EA Price ₩     | 877,200          |



| Cat. No. | C05-79-653 |
|----------|------------|
| Model    | BRK 5436   |



| Cat. No. | C05-79-553 |
|----------|------------|
| Model    | BRK 5448   |



| Cat. No. | C05-79-655 |
|----------|------------|
| Model    | BRK 5494   |

\* 자세한 Spec.은 226 Page 참조