

Advanced Medium Centrifuge, Pro-Analytical 고급형 중형 원심분리기, Max. 1 L

Cat. No.	Model	Type	EA Price ₩
C05-79-175	CR4000	Ambient Medium Centrifuge	4,017,000
C05-79-180	CR4000R	Refrigerated Medium Centrifuge	7,960,100

- LED Display
- 10 RPM 단위 속도조절
- 60 db 이하의 소음
- Emergency Lid-Release
- Rotor 별매
- Max. 1L 까지 가능 (Rotor BRK1000, Bucket B5250 사용시)
- Centurion Scientific Ltd.
- 10 Memory Program
- 0 ~ 99 min 타이머 (30초 단위 조절가능)
- Multi-Point Lid-Locking
- CE / ISO9001 / IVD Certified
- -9 ~ 40°C PID Controlled - CR4000R

Cat. No.	C05-79-175	C05-79-180
Model	CR4000	CR4000R
Size (mm)	w400×d500×h310 / 31 kg	w450×d635×h315 / 61 kg
Speed	500 ~ 15,000 rpm	
Max. Rcf	22,000 G	
Power	310 watts	690 watts



- Reducer는 230 Page 참조
- BRK1000 Rotor는 227 Page 참조



Optional Rotor, Autoclavable at 121°C



Cat. No.	C05-79-549
Model	BRK 5424
Rotor Type	24×2 ml
Speed	500 ~ 15,000 rpm
Tube Size Max.	11×50 mm
Max. Rcf	22,000 G
Radius Max	8.5
Tube Angle	45°
EA Price ₩	877,200



Cat. No.	C05-79-653
Model	BRK 5436
Rotor Type	36×0.5 ml
Speed	500 ~ 15,000 rpm
Tube Size Max.	8×30 mm
Max. Rcf	22,000 G
Radius Max	8.5
Tube Angle	45°
EA Price ₩	971,200



Cat. No.	C05-79-553
Model	BRK 5448
Rotor Type	48×0.2 ml
Speed	500 ~ 15,000 rpm
Tube Size Max.	6×40 mm
Max. Rcf	22,000 G
Radius Max	8.5
Tube Angle	45°
EA Price ₩	930,800



Cat. No.	C05-79-655
Model	BRK 5494
Rotor Type	4×PCR Strips
Speed	500 ~ 15,000 rpm
Tube Size Max.	6×40 mm
Max. Rcf	22,000 G
Radius Max	8.5
Tube Angle	45°
EA Price ₩	971,200



Cat. No.	C05-79-639
Model	BRK 5324
Rotor Type	24×15 ml
Speed	500 ~ 6,000 rpm
Tube Size Max.	17×120 mm
Max. Rcf	4,800 G
Radius Max	12
Tube Angle	30°
EA Price ₩	1,274,200



Cat. No.	C05-79-637
Model	BRK 5308
Rotor Type	8×50 ml
Speed	500 ~ 6,000 rpm
Tube Size Max.	30×120 mm
Max. Rcf	4,800 G
Radius Max	12
Tube Angle	30°
EA Price ₩	1,228,100



Cat. No.	C05-79-635
Model	BRK 5100
Rotor Type	6×100 ml
Speed	500 ~ 6,000 rpm
Tube Size Max.	45×125 mm
Max. Rcf	4,800 G
Radius Max	12
Tube Angle	30°
EA Price ₩	1,405,300



Cat. No.	C05-79-665
Model	BRK 5540
Rotor Type	4×Microtiter Plate
Speed	500 ~ 4,000 rpm
Tube Size Max.	85×128 mm
Max. Rcf	2,500 G
Radius Max	14
Tube Angle	0°
EA Price ₩	2,553,400



Cat. No.	C05-79-659
Model	BRK 5508L
Rotor Type	8×15 ml
Speed	500 ~ 4,000 rpm
Tube Size Max.	17×125 mm
Max. Rcf	2,600 G
Radius Max	14.6
Tube Angle	0°
EA Price ₩	990,800



Cat. No.	C05-79-545
Model	BRK 5401
Rotor Type	24 Capillary & 12×2 ml
Speed	500 ~ 12,000 rpm
Tube Size Max.	2×75 / 11×40 mm
Max. Rcf	13,500 G
Radius Max	8.5
Tube Angle	0 & 60°
EA Price ₩	930,800



Cat. No.	C05-79-563
Model	BRK 5212
Rotor Type	12×15 ml
Speed	500 ~ 10,000 rpm
Tube Size Max.	17×120 mm
Max. Rcf	10,600 G
Radius Max	9.5
Tube Angle	30°
EA Price ₩	2,088,700



Cat. No.	C05-79-565
Model	BRK 5206
Rotor Type	6×50 ml
Speed	500 ~ 10,000 rpm
Tube Size Max.	30×120 mm
Max. Rcf	10,600 G
Radius Max	9.5
Tube Angle	30°
EA Price ₩	1,856,000