

Advanced Large Centrifuge, Pro-Analytical 고급형 대형 원심분리기, Max. 3 L

Cat. No.	Model	Type	EA Price ₩
C05-79-185	CR7000	Ambient Large Centrifuge	9,062,800
C05-79-190	CR7000R	Refrigerated Large Centrifuge	14,092,600

- LED Display
- 10 RPM 단위 속도조절
- 60 db 이하의 소음
- Emergency Lid-Release
- Rotor 별매
- Max. 3L 까지 가능 (Rotor BRK3000, Bucket B3000 사용시)
- Centurion Scientific Ltd.
- 10 Memory Program
- 0 ~ 99 min 타이머 (30초 단위 조절가능)
- Multi-Point Lid-Locking
- CE / ISO9001 / IVD Certified
- -9 ~ 40°C PID Controlled - CR7000R

Cat. No.	C05-79-185	C05-79-190
Model	CR7000	CR7000R
Size (mm)	w600×d630×h375/ 63.5 kg	w998×d630×h410 / 110 kg
Speed	500 ~ 15,000 rpm	
Max. Rcf	22,000 G	
Power	750 watts	1200 watts



- Reducer는 230 Page 참조
- BRK5424 Rotor는 228 Page 참조

Optional Rotor, Autoclavable at 121°C



Cat. No.	C05-79-667
Model	BRK 3020
Rotor Type	4×500 ml
Speed	500 ~ 4,000 rpm
B2000 Bucket 포함	
Max. Rcf	3,600 G
Radius Max	20.5
Tube Angle	0°
EA Price ₩	6,007,600



Cat. No.	C05-79-665
Model	BRK 5540
Rotor Type	4× Microtiter Plate
Speed	500 ~ 4,000 rpm
Tube Size Max.	85×128 mm
Max. Rcf	2,500 G
Radius Max	14
Tube Angle	0°
EA Price ₩	2,553,400



Cat. No.	C05-79-639
Model	BRK 5324
Rotor Type	24×15 ml
Speed	500 ~ 6,000 rpm
Tube Size Max.	17×120 mm
Max. Rcf	4,800 G
Radius Max	12
Tube Angle	30°
EA Price ₩	1,274,200



Cat. No.	C05-79-727
Model	BRK 5224
Rotor Type	24×15 ml
Speed	500 ~ 10,000 rpm
Tube Size Max.	17×120 mm
Max. Rcf	13,400 G
Radius Max	12
Tube Angle	30°
EA Price ₩	2,553,900



Cat. No.	C05-79-725
Model	BRK 5208
Rotor Type	8×50 ml
Speed	500 ~ 10,000 rpm
Tube Size Max.	30×120 mm
Max. Rcf	13,400 G
Radius Max	12
Tube Angle	30°
EA Price ₩	2,553,900



Cat. No.	C05-79-723
Model	BRK 5210
Rotor Type	6×100 ml
Speed	500 ~ 10,000 rpm
Tube Size Max.	45×125 mm
Max. Rcf	13,400 G
Radius Max	12
Tube Angle	30°
EA Price ₩	2,809,300



Cat. No.	C05-79-545
Model	BRK 5401
Rotor Type	24 Capillary & 12×2 ml
Speed	500 ~ 12,000 rpm
Tube Size Max.	2×75 / 11×40 mm
Max. Rcf	13,500 G
Radius Max	8.5
Tube Angle	0 & 60°
EA Price ₩	930,800



Cat. No.	C05-79-653
Model	BRK 5436
Rotor Type	36×0.5 ml
Speed	500 ~ 15,000 rpm
Tube Size Max.	8×30 mm
Max. Rcf	22,000 G
Radius Max	8.5
Tube Angle	45°
EA Price ₩	971,200



Cat. No.	C05-79-553
Model	BRK 5448
Rotor Type	48×0.2 ml
Speed	500 ~ 15,000 rpm
Tube Size Max.	6×40 mm
Max. Rcf	22,000 G
Radius Max	8.5
Tube Angle	45°
EA Price ₩	930,800



Cat. No.	C05-79-655
Model	BRK 5494
Rotor Type	4×PCR Strips
Speed	500 ~ 15,000 rpm
Tube Size Max.	6×40 mm
Max. Rcf	22,000 G
Radius Max	8.5
Tube Angle	45°
EA Price ₩	971,200



Cat. No.	C05-79-637
Model	BRK 5308
Rotor Type	8×50 ml
Speed	500 ~ 6,000 rpm
Tube Size Max.	30×120 mm
Max. Rcf	4,800 G
Radius Max	12
Tube Angle	30°
EA Price ₩	1,228,100



Cat. No.	C05-79-635
Model	BRK 5100
Rotor Type	6×100 ml
Speed	500 ~ 6,000 rpm
Tube Size Max.	45×125 mm
Max. Rcf	4,800 G
Radius Max	12
Tube Angle	30°
EA Price ₩	1,405,300