

Cold Lab Chamber / Refrigerator 저온 냉장 챔버 / 정밀형 냉장고

배송비 별도

Cat. No.	Model	Capacity	Temp. Range	EA Price ₩
C07-02-030	LR-CL185	700 L	0 ~ 10℃	6,297,100
C07-02-040	LR-CL212	1200 L		7,821,400

- 0.01℃ 정밀 제어 구현
- 히터를 통한 정밀한 온도 제어
- ICH 의약품 안정성 시험 가이드라인 충족 (Q1A(R2))
- 내부 순환 Fan을 통한 정밀한 온도 분포도
- Door Open시 Fan 정지 기능, 온도 변화 최소화
- Pair Glass Door, 성에 발생 억제
- LED Lamp를 장착하여 내부관찰 용이
- Graphic LCD Display, 간편하고 직관적인 Controller 조작
- Auto-tuning, Temperature Calibration 기능
- 레코더를 장착하여 온도 체크 가능 (Option)
- 온도편차 알람 기능
- 정전 보상 기능
- Wifi 온도 모니터링 시스템 (Option)
 - 온도 이탈 발생 시 스마트폰으로 알람 발송
 - 설치 장소에 무선 공유기 (wifi) 필요, 스마트폰에 앱 설치 필요
- Validation Service (IQ, OQ, PQ)
- **PL Insurance (제조물 배상 책임 보험)**
- LK Labkorea

Validation Service

Cat. No.	Model	Description	EA Price ₩
C07-91-999	VA-CLC100	Validation for Cold Lab Chamber	

- 당사 기본 양식에 따른 IQ/OQ Validation Service 비용이며, 고객 요청에 의하여 비용이 변경 될 수 있음

Optional Accessory

Cat. No.	Model	Description	EA Price ₩
C07-02-510	TR2002	Compact Paper Recorder	783,900
C07-91-090	RMAL01	Auto Remote Alarm	551,300
C07-02-044	LR-CLSH01	Shelf for LR-CL185	79,000
C07-02-047	LR-CLSH02	Shelf for LR-CL212	88,800

Cat. No.	C07-02-030	C07-02-040	
Model	LR-CL185	LR-CL212.	
Capacity	700 L	1200 L	
Controller	Control	Digital Control	
	Display	GLCD (Graphic LCD)	
	Resolution	0.01℃	
Temperature	Range	0 to +10℃	
	Accuracy at 4℃	±0.3℃	
	Uniformity at 4℃	±1.45℃	±1.89℃
Size (w×d×h)	Interior (mm)	700×700×1400	1200×700×1400
	Exterior (mm)	860×1000×1895	1350×1000×1895
Electric Supply	Power	1 Phase, 220 Vac, 60 Hz	
	Max Consumption	1 kW (4.6 A)	1.5 kW (6.8 A)
	Power Line	Standard Plug	
Material	Interior	Stainless Steel 304	
	Exterior	Powder Coated Steel	
Other	Door	One Door	Two Door
	Wire Shelf	4ea	8ea
	Lamp	LED Lamp × 1ea	
	Caster	Leveling Caster	
Safety	Current Protection, Automatic Fault Sensing System		




Controller



내부이미지



Paper Recorder



Epoxy Coated Wire Shelf