

## Constant Temperature Chamber

### 항온 챔버

배송비 별도

Cat. No.	Model	Capacity	Temp. Range	EA Price ₩
C07-02-300	LI-CTC154P	150 L	-40 ~ +180°C	13,094,400
C07-02-315	LI-CTC254P	250 L		15,778,600
C07-02-320	LI-CTC494P	490 L		18,114,000

- Program Touch Screen Controller
- 가열 속도와 냉각 속도가 빠름. 온도 환경 테스트에 적합
- Max 120 Pattern / 1200 Segment 프로그램 설정 가능
- 운전 중 온도 그래프 표시 기능
- SD메모리 카드를 통한 외부 메모리 지원
- Cable Port 장착, 편리한 테스트 환경 제공
- 내부관찰용 Pair Glass와 LED Lamp 장착
- Door Open시 Fan 정지 기능, 온도 변화 최소화 및 사용자의 안정성 고려
- 레코더를 장착하여 온도 체크 가능 (Option)
- 온도편차 알람 기능
- 정전 보상기능
- Auto-tuning, Temperature Calibration 기능
- 2중화 과온 안전장치, 1차 Controller 과온 Alarm과 2차 과온 차단 회로
- Sign Tower Lamp 장착 가능 (Option)
- Industrial Plug 적용 제품, 구매 전 확인 필요
- Validation Service (IQ, OQ, PQ)
- **PL Insurance (제조물 배상 책임 보험)**
- LK Labkorea

### Optional Accessory

Cat. No.	Model	Description	EA Price ₩
C07-02-500	TRL-6	Paper Recorder	2,223,100
C07-02-341	LI-CTCSH01	Shelf for 150 L (LI-CTC154P)	103,600
C07-02-343	LI-CTCSH02	Shelf for 250 L (LI-CTC254P)	111,900
C07-02-344	LI-CTCSH04	Shelf for 490 L (LI-CTC494P)	120,100
O01-96-100	SI300	Sign Tower Lamp	271,500

### Validation Service

Cat. No.	Model	Description	EA Price ₩
C07-91-977	VA-TC100	Validation for Constant Temp. Chamber	

- 당사 기본 양식에 따른 IQ/OQ Validation Service 비용이며, 지역이나 고객 요청사항에 따라 비용이 변경될 수 있음



PL



Program Temp. Controller



Cable Port




Plate Shelf



Paper Recorder



SI300

Cat. No.	C07-02-300	C07-02-315	C07-02-320	
Model	LI-CTC154P	LI-CTC254P	LI-CTC494P	
Capacity	150 L	250 L	490 L	
Controller	Control	Program Touch Screen Control		
	Display	5.7" Color LCD		
	Resolution	0.01°C		
Temperature	Range	-40 ~ +180°C		
	Accuracy at 50°C	±0.3°C		
	Uniformity at 50°C	±0.82°C	±1.13°C	±1.88°C
Size	Interior	w500 × d500 × h600 mm	w600 × d600 × h700 mm	w700 × d700 × h1000 mm
	Exterior	w720 × d940 × h1730 mm	w820 × d1040 × h1830 mm	w920 × d1140 × h2200 mm
Material	Interior	Stainless Steel 304		
	Exterior	Powder Coated Steel		
Electric Supply	Power	1 Phase / 220VAC / 60 Hz		3 Phase / 380VAC / 60 Hz
	Max Consumption	3 kW (13.6 A)	4.3 kW (19.5 A)	7.3 kW (11.1 A)
	Power Line	Industrial Plug		
Other	Wire Shelf	SUS304 Punching Plate Shelves / 2ea		3 ea
	Cable Port	Ø50 mm × 1ea		
	Caster	Leveling Caster		
Safety	Over Temp. & Current Protection, Emergency switch			