

Chromatography Column 크로마토그래피 컬럼, PTFE 콕크, LukeGL®

- Borosilicate Glass 3.3
- PTFE-Plug Stopcock, Straight Bore
- Without Filter Disc
- 요청에 의해 규격 외 사이즈도 제작 공급 가능함
- LukeGL®

Column, 기본형

Cat. No.	Model	Bore	용량	컬럼 od	길이 (L)	EA Price ₩	
C08-10-006	5435701	2.5 mm	7 ml	Ø10 mm	200 mm	38,000	
C08-10-012	5435702		10 ml		250 mm	39,200	
C08-10-018	5435703		12 ml		300 mm	39,200	
C08-10-024	5435704		23 ml	Ø15 mm	250 mm	41,000	
C08-10-030	5435705		28 ml		300 mm	42,200	
C08-10-036	5435706		38 ml	Ø19 mm	400 mm	42,200	
C08-10-042	5435707		44 ml		250 mm	41,000	
C08-10-048	5435708		52 ml		300 mm	42,200	
C08-10-054	5435709		70 ml	Ø25 mm	400 mm	42,200	
C08-10-060	5435710		85 ml		250 mm	44,500	
C08-10-066	5435711		100 ml	Ø38 mm	300 mm	44,500	
C08-10-072	5435712		135 ml		400 mm	45,400	
C08-10-078	5435713		170 ml		500 mm	45,400	
C08-10-084	5435714		270 ml	Ø50 mm	300 mm	51,800	
C08-10-090	5435715		360 ml		400 mm	55,200	
C08-10-096	5435716		450 ml		500 mm	58,700	
C08-10-102	5435717		4.0 mm	650 ml	Ø65 mm	400 mm	65,100
C08-10-108	5435718			830 ml		500 mm	68,700
C08-10-114	5435719	990 ml		600 mm		70,400	
C08-10-120	5435720	6.0 mm	1130 ml	Ø75 mm	400 mm	77,100	
C08-10-126	5435721		1690 ml		600 mm	83,100	
C08-10-132	5435722		1530 ml	Ø50 mm	400 mm	94,000	
C08-10-138	5435723		2300 ml		600 mm	114,500	

Column with Top Joint, 상부 조인트 부착형 (Joint Size 24/40)

Cat. No.	Model	Bore	용량	컬럼 od	길이 (L)	EA Price ₩	
C08-10-282	5435901	2.5 mm	7 ml	Ø10 mm	200 mm	47,000	
C08-10-288	5435902		10 ml		250 mm	47,000	
C08-10-294	5435903		12 ml		300 mm	50,600	
C08-10-300	5435904		23 ml	Ø15 mm	250 mm	48,600	
C08-10-306	5435905		28 ml		300 mm	48,800	
C08-10-312	5435906		38 ml	Ø19 mm	400 mm	52,100	
C08-10-318	5435907		44 ml		250 mm	50,100	
C08-10-324	5435908		52 ml		300 mm	50,100	
C08-10-330	5435909		70 ml	Ø25 mm	400 mm	53,600	
C08-10-336	5435910		85 ml		250 mm	51,200	
C08-10-342	5435911		100 ml	Ø38 mm	300 mm	51,200	
C08-10-348	5435912		135 ml		400 mm	53,600	
C08-10-354	5435913		170 ml		500 mm	54,100	
C08-10-360	5435914		270 ml	Ø50 mm	300 mm	61,900	
C08-10-366	5435915		360 ml		400 mm	62,200	
C08-10-372	5435916		450 ml		500 mm	65,300	
C08-10-378	5435917		4.0 mm	650 ml	Ø65 mm	400 mm	73,100
C08-10-384	5435918			830 ml		500 mm	84,300
C08-10-390	5435919	990 ml		600 mm		85,800	
C08-10-396	5435920	6.0 mm	1130 ml	Ø75 mm	400 mm	101,200	
C08-10-402	5435921		1690 ml		600 mm	106,900	
C08-10-408	5435922		1530 ml	Ø50 mm	400 mm	138,600	
C08-10-414	5435923		2300 ml		600 mm	156,600	

- Joint Size는 요청에 의해 변경 가능함

