

## Chromatography Column with Glass Filter 크로마토그래피 컬럼, PTFE 콕크, LukeGL®

- Borosilicate Glass 3.3
- PTFE-Plug Stopcock, Straight Bore
- Filter Pore Size : 40 ~ 60 um (요청에 의해 변경 가능)
- 요청에 의해 규격 외 사이즈도 제작 공급 가능함
- LukeGL®

### Column, 기본형

Cat. No.	Model	Bore	용량	컬럼 od	길이 (L)	EA Price ₩
C08-10-144	5435601	2.5 mm	7 ml	Ø10 mm	200 mm	47,000
C08-10-150	5435602		10 ml		250 mm	49,000
C08-10-156	5435603		12 ml		300 mm	49,000
C08-10-162	5435604		23 ml	Ø15 mm	250 mm	50,100
C08-10-168	5435605		28 ml		300 mm	50,100
C08-10-174	5435606		38 ml		400 mm	50,600
C08-10-180	5435607		44 ml	Ø19 mm	250 mm	50,100
C08-10-186	5435608		52 ml		300 mm	53,400
C08-10-192	5435609		70 ml		400 mm	53,400
C08-10-198	5435610		85 ml	Ø25 mm	250 mm	62,000
C08-10-204	5435611		100 ml		300 mm	62,000
C08-10-210	5435612		135 ml		400 mm	66,400
C08-10-216	5435613	170 ml	500 mm		66,400	
C08-10-222	5435614	270 ml	Ø38 mm		300 mm	77,700
C08-10-228	5435615	360 ml			400 mm	77,700
C08-10-234	5435616	450 ml	Ø50 mm	500 mm	81,000	
C08-10-240	5435617	650 ml		400 mm	90,000	
C08-10-246	5435618	830 ml		500 mm	102,400	
C08-10-252	5435619	990 ml	Ø65 mm	600 mm	108,400	
C08-10-258	5435620	1130 ml		400 mm	118,000	
C08-10-264	5435621	1690 ml		600 mm	134,900	
C08-10-270	5435622	6.0 mm	1530 ml	Ø75 mm	400 mm	151,800
C08-10-276	5435623		2300 ml		600 mm	166,000

### Column with Top Joint, 상부 조인트 부착형 (Joint Size 24/40)

Cat. No.	Model	Bore	용량	컬럼 od	길이 (L)	EA Price ₩
C08-10-420	5435801	2.5 mm	7 ml	Ø10 mm	200 mm	55,400
C08-10-426	5435802		10 ml		250 mm	57,300
C08-10-432	5435803		12 ml		300 mm	58,600
C08-10-438	5435804		23 ml	Ø15 mm	250 mm	59,400
C08-10-444	5435805		28 ml		300 mm	59,600
C08-10-450	5435806		38 ml		400 mm	59,600
C08-10-456	5435807		44 ml	Ø19 mm	250 mm	63,500
C08-10-462	5435808		52 ml		300 mm	64,100
C08-10-468	5435809		70 ml		400 mm	64,100
C08-10-474	5435810		85 ml	Ø25 mm	250 mm	68,300
C08-10-480	5435811		100 ml		300 mm	68,300
C08-10-486	5435812		135 ml		400 mm	72,000
C08-10-492	5435813	170 ml	500 mm		72,500	
C08-10-498	5435814	270 ml	Ø38 mm		300 mm	80,400
C08-10-504	5435815	360 ml			400 mm	84,300
C08-10-510	5435816	450 ml	Ø50 mm	500 mm	87,700	
C08-10-516	5435817	650 ml		400 mm	106,500	
C08-10-522	5435818	830 ml		500 mm	111,000	
C08-10-528	5435819	990 ml	Ø65 mm	600 mm	115,800	
C08-10-534	5435820	1130 ml		400 mm	134,900	
C08-10-540	5435821	1690 ml		Ø75 mm	600 mm	179,400
C08-10-546	5435822	1530 ml	400 mm		187,100	
C08-10-552	5435823	2300 ml	600 mm	197,300		

- Joint Size는 요청에 의해 변경 가능함

