

TLC Plate TLC 플레이트

• UNIPLATE® TLC Plate

- ① Silicagel G : No Organic Substance
 - Water를 20% 이상 포함한 Solution과의 사용은 불가함
 - Compatible with Sulfuric Acid Charring Visualization
 - ② Silicagel GF : "G" + Fluorescent Indicator : 250 nm (UV Short Wave)
- Plate Size : w200×d200 mm • Analtech / UNIPLATE®

Glass Plate, UNIPLATE®

Cat. No.	Model	Baking	형광 처리	CS Price ₩
C08-148-110	01011	Silicagel G	No	25ea 267,200
C08-148-113	02011	Silicagel GF	Yes	25ea 267,200

- 250 um Thickness, 10 um Particle, 60 Å pore

Glass Plate, Preparative, UNIPLATE®

Cat. No.	Model	Baking	두께	형광 처리	CS Price ₩
C08-148-003	01012	Silicagel G	500 um	No	25ea 430,500
C08-148-006	02012	Silicagel GF		Yes	25ea 430,500
C08-148-009	01013	Silicagel G	1000 um	No	25ea 517,800
C08-148-012	02013	Silicagel GF		Yes	25ea 517,800
C08-148-015	01014	Silicagel G	1500 um	No	25ea 654,800
C08-148-018	02014	Silicagel GF		Yes	25ea 654,800
C08-148-021	01015	Silicagel G	2000 um	No	25ea 770,200
C08-148-024	02015	Silicagel GF		Yes	25ea 770,200

- Large Quantify Material의 Separation용

Aluminum Plate, UNIPLATE®

Cat. No.	Model	Baking	형광 처리	CS Price ₩
C08-148-506	156017	Silicagel G	No	25ea 340,100
C08-148-509	157017	Silicagel GF	Yes	25ea 340,100

- 200 um Thickness

Polyester Plate, UNIPLATE®

Cat. No.	Model	Baking	형광 처리	CS Price ₩
C08-148-706	158017	Silicagel G	No	25ea 340,100
C08-148-709	159017	Silicagel GF	Yes	25ea 340,100

- 200 um Thickness




Chromatography Paper / Blotting Paper 크로마토 그래피 페이퍼

Cat. No.	Model	Grade	두께	PK Price ₩
C08-07-050	3030917	3 MM Chr	0.34 mm	100ea 488,000
C08-07-073	3017915	17 Chr	0.92 mm	250ea 571,000
C08-07-100	3001917	1 Chr	0.18 mm	100ea 383,000

- 재질 : Pure Cellulose • Size : w460×d570 mm • Whatman

용도

Grade	유속	용도
3 MM Chr	130 mm/30min	Extensively for General Chromatography and Electrophoresis
17 Chr	190 mm/30min	The Heaviest Loadings and for Preparative Paper Chromatography and Electrophoresis
1 Chr	130 mm/30min	General Analytical Separation



TLC Spotting Guide TLC 스포팅 가이드

Cat. No.	Model	Description	EA Price ₩
C08-148-901	25-00	up to 20×20 cm Plate	88,000

- with Metric Scale
- Analtech

