


## Mini TLC Developing Chamber 소형 TLC 전개조

Cat. No.	Model	Description	Size	EA Price ₩
C08-132-007	1184-61	10 cm Plate용	w120 x d64 x h115mm	131,600
C08-132-010	1184-62	TLC 전개조	w120 x d86 x h115mm	153,500

- 소형 유리 TLC 전개조, 뚜껑 포함
- 10×10 cm TLC Plate용
- 내부 폭 86mm인 1184-62는 TLC Plate Rack(1185-10)과 같이 사용 가능
- Glass Blowing

### Optional Accessory

Cat. No.	Model	Description	EA Price ₩
C08-132-014	1184-63	1184-61용 Lid	18,100
C08-132-018	1184-65	1184-62용 Lid	18,100



1184-61

## TLC Developing Chamber TLC 전개조

Cat. No.	Model	Description	EA Price ₩
C08-132-021	1184-01	20 cm Plate용 TLC 전개조	312,000
C08-132-028	1184-03	Replacement Lid Only	39,600

- 20×20 cm TLC Plate용
- 내부 Size : w270×d70×h260 mm
- 뚜껑 포함
- Chemglass



1184-01

## Latch-Lid TLC Developing Chamber 덮개 고정식 TLC 전개조

Cat. No.	Model	Description	EA Price ₩
C08-132-031	1183-01	20 cm Plate용 TLC 전개조	451,000
C08-132-038	1183-03	Replacement Lid Only	92,800

- 20×20 cm TLC Plate용
- 내부 Size : w270×d70×h260 mm
- 고정식 덮개 포함
- TLC Plate Rack은 별도 구매하여야 함
- Chemglass



1183-01

## TLC Plate Rack TLC 전개조 랙

Cat. No.	Model	사용 TLC Chamber Model	TLC Plate	EA Price ₩
C08-132-045	1185-10	1184-62	wide 10cm용	407,600
C08-132-049	1185-20	1183-01, 1184-01	wide 20cm용	499,300

- 재질 : Teflon Coated Aluminum
- TLC Plate 거치용
- 1185-10 : 10x10cm 또는 10x5cm TLC Plate용
- 1185-20 : 20x20cm TLC Plate용
- Capacity : 6장 (3.1t), 또는 12장 (1.59t)
- 1185-10은 1184-62 용으로만 사용 가능
- Chemglass



1185-20

- Kimble 소형 TLC Chamber는 62 Page 참조

