

## SPE Cartridge SPE 카트리지

- GC, HPLC 등을 통해 분석 전 시료 정화 및 농축 가능
- 2개의 20um polyethylene frit이 있는 고순도 polypropylene body
- Multi-position 진공 매니폴드 사용을 위한 male Luer-tip 장착
- 4단계 Sample preparation (Conditioning, Sample Application, Wash, Elution) 가능
- JGF (J. G. Finneran)

### 1ml SPE Cartridge, 100mg Dosage

Cat. No.	Model	Packing	Miligrams	CS Price ₩
C08-214-011	JGF01C1801	C18	100	100ea 368,600
C08-214-015	JGF01C801	C8		100ea 386,300
C08-214-019	JGF01CN01	CN		100ea 388,900
C08-214-023	JGF01SI01	Silica		100ea 311,200
C08-214-027	JGF01NH01	NH2		100ea 388,900
C08-214-031	JGF01FL0160	Florisil® (60/100)		100ea 317,200
C08-214-035	JGF01FL01100	Florisil® (100/200)		100ea 317,200
C08-214-039	JGF01DI01	Diol		100ea 388,900
C08-214-043	JGF01SA01	SAX		100ea 452,000
C08-214-047	JGF01SC01	SCX		100ea 465,900

### 3ml SPE Cartridge, 500mg Dosage

Cat. No.	Model	Packing	Miligrams	CS Price ₩
C08-214-101	JGF03C1805	C18	500	50ea 344,000
C08-214-105	JGF03C805	C8		50ea 411,700
C08-214-109	JGF03CN05	CN		50ea 392,000
C08-214-113	JGF03SI05	Silica		50ea 204,600
C08-214-117	JGF03NH05	NH2		50ea 392,000
C08-214-121	JGF03FL0560	Florisil® (60/100)		50ea 215,500
C08-214-125	JGF03FL05100	Florisil® (100/200)		50ea 215,500
C08-214-129	JGF03DI05	Diol		50ea 392,000
C08-214-133	JGF03SA05	SAX		50ea 551,700
C08-214-137	JGF03SC05	SCX		50ea 586,500

### 6ml SPE Cartridge, 1000mg Dosage

Cat. No.	Model	Packing	Miligrams	CS Price ₩
C08-214-191	JGF06C1810	C18	1000	30ea 331,800
C08-214-195	JGF06C810	C8		30ea 411,200
C08-214-199	JGF06CN10	CN		30ea 387,000
C08-214-203	JGF06SI10	Silica		30ea 160,700
C08-214-207	JGF06NH10	NH2		30ea 387,000
C08-214-211	JGF06FL1060	Florisil® (60/100)		30ea 182,700
C08-214-215	JGF06FL10100	Florisil® (100/200)		30ea 182,700
C08-214-219	JGF06DI10	Diol		30ea 386,300
C08-214-223	JGF06SA10	SAX		30ea 562,500
C08-214-227	JGF06SC10	SCX		30ea 620,700

## SPE Cartridge SPE 카트리지

Cat. No.	Model	Packing	Miligrams	PK Price ₩
C08-148-822	01-00	Silica Gel	300	48ea 444,200
C08-148-825	01-08	UNIBOND™ C2	400	48ea 600,400
C08-148-828	01-09	UNIBOND™ C8		48ea 600,400
C08-148-831	01-10	UNIBOND™ C18		48ea 600,400
C08-148-834	01-22	UNIBOND™ Amino	1000	48ea 600,400
C08-148-837	01-23	UNIBOND™ Cyano		48ea 600,400
C08-148-840	01-33	Alumina (Neutral)		48ea 541,000
C08-148-843	01-34	Alumina (Basic)	500	48ea 541,000
C08-148-846	01-35	Alumina (Acidic)		48ea 541,000
C08-148-849	01-96	Florisil™	500	48ea 479,400

- Prepacked with Adsorbent and Ready for Use
- 구성 : Medical Grade PP Cartridge, 1 cc Packing
- Porous Frit Outlet Filter Combination
- Remove Sample Particulate down to 10 um



### 3ml SPE Cartridge, 200mg Dosage

Cat. No.	Model	Packing	Miligrams	CS Price ₩
C08-214-055	JGF03C1802	C18	200	50ea 233,600
C08-214-059	JGF03C802	C8		50ea 252,500
C08-214-063	JGF03CN02	CN		50ea 248,100
C08-214-067	JGF03SI02	Silica		50ea 228,900
C08-214-071	JGF03NH02	NH2		50ea 245,100
C08-214-075	JGF03FL0260	Florisil® (60/100)		50ea 172,500
C08-214-079	JGF03FL02100	Florisil® (100/200)		50ea 172,500
C08-214-083	JGF03DI02	Diol		50ea 245,100
C08-214-087	JGF03SA02	SAX		50ea 300,200
C08-214-091	JGF03SC02	SCX		50ea 324,200

### 6ml SPE Cartridge, 500mg Dosage

Cat. No.	Model	Packing	Miligrams	CS Price ₩
C08-214-145	JGF06C1805	C18	500	30ea 210,100
C08-214-149	JGF06C805	C8		30ea 250,600
C08-214-153	JGF06CN05	CN		30ea 238,800
C08-214-157	JGF06SI05	Silica		30ea 126,400
C08-214-161	JGF06NH05	NH2		30ea 238,800
C08-214-165	JGF06FL0560	Florisil® (60/100)		30ea 130,600
C08-214-169	JGF06FL05100	Florisil® (100/200)		30ea 130,600
C08-214-173	JGF06DI05	Diol		30ea 238,800
C08-214-177	JGF06SA05	SAX		30ea 320,400
C08-214-181	JGF06SC05	SCX		30ea 341,600

### 12ml SPE Cartridge, 1000mg Dosage

Cat. No.	Model	Packing	Miligrams	CS Price ₩
C08-214-235	JGF12C1810	C18	1000	20ea 258,400
C08-214-239	JGF12C810	C8		20ea 289,400
C08-214-243	JGF12CN10	CN		20ea 296,700
C08-214-247	JGF12SI10	Silica		20ea 151,700
C08-214-251	JGF12NH10	NH2		20ea 296,700
C08-214-255	JGF12FL1060	Florisil® (60/100)		20ea 165,400
C08-214-259	JGF12FL10100	Florisil® (100/200)		20ea 165,400
C08-214-263	JGF12DI10	Diol		20ea 296,700
C08-214-267	JGF12SA10	SAX		20ea 420,700
C08-214-271	JGF12SC10	SCX		20ea 465,900



- Connection : (Inlet) Female Lock, (Outlet) Male Slip
- Syringe 및 Vacuum Manifold와 같이 사용가능
- Typical Flow Rate : 10 cc/min at 14 psi
- Analtech / SPICE®