

2-Prong Extension Clamp / 2-P Clamp 2-P 클램프

Cat. No.	Model	Max. Grip	길이	EA Price ₩
C09-143-005	L2P-40	60 mm	250 mm	7,800
C09-143-006	L2P-40L		605 mm	11,900

- 다이캐스팅후 크롬 도금 • Rod : Ø12.7 mm
- 미끄럼 방지를 위한 그립부 비닐코팅
- LK Labkorea



2-Prong Extension Clamp / 2-P Clamp 2-P 클램프

Cat. No.	Model	Grip	길이	EA Price ₩
C09-284-001	2544	12 ~ 25 mm	235 mm	25,700
C09-284-004	2554	20 ~ 40 mm	240 mm	27,500
C09-284-007	2564	38 ~ 60 mm	253 mm	29,000
C09-284-010	2574	30 ~ 80 mm	257 mm	32,300


- Nickel Plated, Vinyl Coated • Rod : Ø10 mm
- Usbeck



3-Prong Extension Clamp / 3-P Clamp 3-P 클램프

Cat. No.	Model	규격	Max. Grip	길이	EA Price ₩
C09-143-010	L3P-20	소	30 mm	210 mm	7,200
C09-143-015	L3P-40	중	60 mm	245 mm	7,700
C09-143-020	L3P-70	대	80 mm	270 mm	8,300
C09-143-022	L3P-70L	대/장형		615 mm	11,900


- 다이캐스팅후 크롬 도금 • 미끄럼 방지를 위한 그립부 비닐코팅
- 2-P 클램프보다 견고하고 안정적인
- LK Labkorea



3-Prong Extension Clamp / 3-P Clamp 3-P 클램프

Cat. No.	Model	Grip	길이	EA Price ₩
C09-284-013	2801	0 ~ 40 mm	227 mm	17,200
C09-284-017	2802	25 ~ 80 mm	245 mm	27,100

- Zinc Diecasting, Natural Alloy
- Vinyl Coated, Rod : Ø10 mm
- Usbeck



Condenser Clamp 콘덴서 클램프

Cat. No.	Model	규격	Max. Grip	EA Price ₩
C09-503-071	33430302	4-Prong	up to 125mm	29,500

- Aluminium, Powder Coated with Steel Rod
- Silicon Sleeve
- LK Labkorea

