

Boss Head Clamp 보스 헤드 클램프

Cat. No.	Model	Max. Grip	길이	EA Price ₩
H03-06-193	HSRC20	12.8 mm	77 ~ 95 mm	8,700


- Aluminum Anodizing
- 사용 가능한 봉의 직경은 12.8mm 임
- LK Labkorea



Square Holder 사각 홀더

Cat. No.	Model	Pipe Size	재질	EA Price ₩
C09-143-100	LCP-FR01	12.7 mm	크롬 도금 ①	2,900
C09-143-103	LCP-FR11		ABS ②	2,900


- 90° 용 홀더
- 사용 가능한 봉의 직경은 12.7 mm 임
- LK Labkorea



Hook Connector 90° 훅 컨넥터

Cat. No.	Model	Pipe Size	재질	EA Price ₩
C09-143-206	LCP-FR04	12.7 mm	크롬 도금 ①	3,400
C09-143-208	LCP-FR14		Stainless Steel ②	2,700

- LCP-FR04는 Guide Holder가 내장 되어 있음
- Frame 조립시 미끄러지거나 비틀림 없이 안전함
- LK Labkorea



Bosshead Clamp 원형 보스헤드 클램프

Cat. No.	Model	Crossbore	Plate Slot	EA Price ₩
C09-284-020	2424	Ø13 mm	Ø6 mm	42,500

- Aluminum, Nickel Plated Brass
- Pipe, Rod, Plate 고정용으로 사용
- Bores on the face : Ø10mm and Ø13mm
- Usbeck



Bosshead Clamp 보스헤드 클램프

Cat. No.	Model	Rod	Screw	EA Price ₩
C09-284-064	2420	up to Ø13 mm	Thumb-type	25,100
C09-284-067	2421		Hinged-type	30,500

- Aluminum, Nickel Plated Brass
- Nickel Plated Brass Screw
- Usbeck



Globular Bosshead Clamp 구형 보스헤드 클램프

Cat. No.	Model	Rod	재질	EA Price ₩
C09-284-070	2360	up to Ø13 mm	Powder Coated	15,400
C09-284-073	2362		Natural zinc diecasting	14,300
C09-284-076	2363		Nickel Plated	18,200

- Zinc Diecasting, Hinged Type Nickel Plated
- Body Size : Ø35 mm
- Usbeck

