

**Plastic Conical Joint Clamp, "KECK®-clips"**  
**플라스틱 코니칼 조인트 클램프**

Cat. No.	Model	규격	색상	EA Price ₩
C09-83-006	KC10	ST10	Green	5,200
C09-83-012	KC12	ST12	Violet	5,300
C09-83-018	KC14	ST14	Yellow	5,500
C09-83-024	KC19	ST19	Blue	5,600
C09-83-030	KC24	ST24	Green	6,000
C09-83-036	KC29	ST29	Red	6,400
C09-83-042	KC34	ST34	Orange	7,000
C09-83-048	KC40	ST40	Yellow	8,000
C09-83-054	KC45	ST45	Brown	8,700

- 재질 : POM, -40 ~ 130°C 내열      • 강한 텐션과 안정성 및 내화학성이 있음
- "Made in Swiss" 의 고급형 제품
- DWK Lifescience



**Plastic Conical Joint Clamp**  
**경제형 플라스틱 조인트 클램프**

Cat. No.	Model	규격	색상	EA Price ₩
C09-339-006	KC10-DO	ST10	Green	2,050
C09-339-012	KC12-DO	ST12	Violet	2,050
C09-339-018	KC14-DO	ST14	Yellow	2,170
C09-339-024	KC19-DO	ST19	Blue	2,530
C09-339-030	KC24-DO	ST24	Green	2,530
C09-339-036	KC29-DO	ST29	Red	2,650
C09-339-042	KC34-DO	ST34	Orange	3,250
C09-339-048	KC40-DO	ST40	Yellow	3,860
C09-339-054	KC45-DO	ST45	Brown	3,980

- 재질 : POM, -40 ~ 130°C 내열      • 강한 텐션과 안정성 및 내화학성이 있음
- 보급형 코니칼 조인트 클램프
- LabPlasti®



**Metal Conical Joint Clamp, "KECK®-clips"**  
**메탈 코니칼 조인트 클램프**

Cat. No.	Model	규격	EA Price ₩	PK Price ₩
C09-83-060	KCM14	ST14	11,700	10ea 109,200
C09-83-066	KCM19	ST19	12,300	10ea 114,300
C09-83-072	KCM24	ST24	12,900	10ea 120,500
C09-83-078	KCM29	ST29	14,200	10ea 132,500

- 재질 : Nikel, 200°C 내열      • POM 플라스틱 보다 강한 텐션과 내열성을 가짐
- DWK Life Science



**Stainless Steel Conical Joint Clamp**  
**스테인레스 코니칼 조인트 클램프**

Cat. No.	Model	규격	EA Price ₩
C09-132-072	KCS24	ST24	8,000

- 재질 : Stainless Steel      • 고온 작업에 알맞은 메탈 조인트 클램프
- LK Labkorea



**Flask Joint Clamp**  
**플라스크 조인트 클램프**

Cat. No.	Model	규격	용량	EA Price ₩
C09-284-079	5170	NS 14	100 ml	5,500
C09-284-082	5171	NS 19	250 ml	6,100
C09-284-085	5172	NS 24	500 ml	6,500
C09-284-088	5173	NS 29	1 L	7,300

- Nickel Plated Steel      • Usbeck

