


Clean Bench, PCR Prep Station
다용도 클린벤치

Cat. No.	Model	Dimension	EA Price ₩
C10-193-401	MY-PCR24-010	w610×d610×h762 mm	6,820,400
C10-193-411	MY-PCR32-010	w813×d610×h762 mm	7,316,300
C10-193-421	MY-PCR48-010	w1219×d610×h762 mm	10,166,000

- PCR 실험 시 샘플의 오염 위험 감소
- Vertical Laminar Flow, HEPA Filter (Class 100)
- Polypropylene Work Surface (화학물 불침투성 우수)
- 피펫 홀더와 통합된 UV 선반
- UV Light 타이머
- 생물안전작업대로 사용 불가
- Mystaire

Spare Accessory

Cat. No.	Model	Description	EA Price ₩
C10-193-441	MY-9001	w24" × d24" × h34" Cart for MY-PCR24	2,408,700
C10-193-445	MY-9002	w32" × d24" × h34" Cart for MY-PCR32	2,729,900
C10-193-449	MY-9004	w48" × d24" × h34" Cart for MY-PCR48	3,051,000
C10-193-453	MY-HEPA18	HEPA filter for MY-PCR24, MY-PCR48	1,043,700
C10-193-457	MY-HEPA	HEPA filter for MY-PCR32	1,124,100
C10-193-461	MY-PRE2	Pre-filter (12/pk) for MY-PCR24	634,200
C10-193-465	MY-PRELF	Pre-filter (12/pk) for MY-PCR32	550,000
C10-193-469	MY-PRE7	Pre-filter (12/pk) for MY-PCR48	1,011,700
C10-193-473	MY-UVB4	Shortwave UV Bulb (4/pk) for MY-PCR24, MY-PCR32	622,300



Cat. No.	C10-193-401	C10-193-411	C10-193-421	
Model	MY-PCR24-010	MY-PCR32-010	MY-PCR48-010	
Size	External	w610×d610×h762 mm	w813×d610×h762 mm	
	Internal	w584×d584×h457 mm	w787×d584×h457 mm	
	Opening (w × Opening1 × Opening2)	533 × 254 × 508 mm	737 × 254 × 508 mm	1143 × 254 × 508 mm
Filter	Main	1ea	1ea	2ea
	Pre-filter	1ea	1ea	2ea
Air Quality	Class 100 / ISO 5			
Blower	Brushless-Sparkless			
Material	Work Surface	Polypropylene		
	Shell	Polycarbonate		
	Front Sash	Polycarbonate		
	Filtration Housing	ABS		
	Pipette Shelf	Polycarbonate		
Electrical	220V, 0.7A, 150W		220V, 1.2A, 270W	
Lighting	Fluorescent			
UV	UV Light	254nm Shortwave		
	UV Light Timer	15 Minutes		