




Graham Condenser 나사관 냉각기, LukeGL®

Cat. No.	Model	길이	Socket / Cone NS	Water Connection	EA Price ₩	
C18-10-006	5240201	150 mm	14/23	Nozzle ①	32,500	
C18-10-008	5240229	200 mm			36,100	
C18-10-012	5240202	150 mm	24/40		35,100	
C18-10-018	5240203	200 mm			36,100	
C18-10-022	5240231	250 mm			41,000	
C18-10-024	5240204	300 mm			41,400	
C18-10-030	5240205	400 mm			46,700	
C18-10-036	5240206	500 mm			50,700	
C18-10-048	5240211	150 mm	14/23		GL14 with PP Connector ②	43,100
C18-10-054	5240212		45,800			
C18-10-060	5240213	200 mm	24/40	47,000		
C18-10-062	5340241	250 mm		50,600		
C18-10-066	5240214	300 mm		51,900		
C18-10-072	5240215	400 mm		58,100		
C18-10-078	5240216	500 mm		62,400		


- Borosilicate Glass 3.3 • Joint Size는 요청에 의해 변경 가능함
- LukeGL®



Dimroth Condenser, Reflux-A, Reservoir 환류 냉각기, LukeGL®

Cat. No.	Model	길이	Socket / Cone NS	Water Connection	EA Price ₩	
C18-10-088	5241200	100 mm	14/23	Nozzle ①	40,000	
C18-10-090	5241201	150 mm			41,200	
C18-10-096	5241202		41,200			
C18-10-102	5241203	200 mm	24/40		41,600	
C18-10-106	5241231	250 mm			44,600	
C18-10-108	5241204	300 mm			44,800	
C18-10-114	5241205	400 mm			51,100	
C18-10-120	5241206	500 mm			57,500	
C18-10-132	5241211	150 mm			14/23	GL14 with PP Connector ②
C18-10-138	5241212		50,800			
C18-10-144	5241213	200 mm	24/40	52,300		
C18-10-146	5241241	250 mm		54,900		
C18-10-150	5241214	300 mm		56,100		
C18-10-156	5241215	400 mm		68,700		
C18-10-162	5241216	500 mm		75,400		

- Borosilicate Glass 3.3 • Joint Size는 요청에 의해 변경 가능함
- LukeGL®



Condenser, Reflux-B 환류 냉각기, LukeGL®

Cat. No.	Model	길이	Socket / Cone NS	Water Connection	EA Price ₩	
C18-10-174	5240301	150 mm	14/23	Nozzle ①	38,600	
C18-10-180	5240302				40,200	
C18-10-186	5240303	200 mm	24/40		47,100	
C18-10-188	5240331	250 mm			49,000	
C18-10-192	5240304	300 mm			50,600	
C18-10-198	5240305	400 mm			54,200	
C18-10-204	5240306	500 mm			59,000	
C18-10-216	5240311	150 mm			14/23	GL14 with PP Connector ②
C18-10-222	5240312		49,900			
C18-10-228	5240313	200 mm	24/40		56,700	
C18-10-230	5240341	250 mm		58,700		
C18-10-234	5240314	300 mm		60,200		
C18-10-240	5240315	400 mm		63,900		
C18-10-246	5240316	500 mm		68,700		

- Borosilicate Glass 3.3 • Joint Size는 요청에 의해 변경 가능함
- LukeGL®

