




### Davis Condenser 이중 자켓 냉각기, LukeGL®

Cat. No.	Model	Jacket 길이	Socket Cone NS	Water Connection	EA Price ₩
C18-10-476	5215303	200 mm	24/40	Nozzle ①	62,400
C18-10-478	5215331	250 mm			82,000
C18-10-480	5215304	300 mm			84,300
C18-10-486	5215305	400 mm			87,700
C18-10-492	52153G3	200 mm		GL14 with PP Connectors ②	79,300
C18-10-494	52153G9	250 mm			100,000
C18-10-498	52153G4	300 mm			102,900
C18-10-504	52153G5	400 mm			104,600

- Borosilicate Glass 3.3
- Joint Size는 요청에 의해 변경 가능함
- LukeGL®



### Dry Ice Condenser 드라이 아이스 냉각기, LukeGL®

- Borosilicate Glass 3.3
- With Cone and Water Connections
- LukeGL®

#### With Hose Nipple ①

Cat. No.	Model	Inside Size		NS Cone	EA Price ₩
		Inner Dia.	Depth		
C18-10-510	L120900	ø30 mm	120 mm	14/23	55,700
C18-10-516	L120901	ø40 mm	200 mm	24/40	68,200
C18-10-522	L120903	ø50 mm	250 mm		84,300
C18-10-528	L120905	ø75 mm	300 mm		134,000

#### With Joint ②

Cat. No.	Model	Inside Size		NS Cone × 2개	EA Price ₩
		Inner Dia.	Depth		
C18-10-534	L1209A22	ø30 mm	120 mm	14/23	64,200
C18-10-540	L1209A26	ø40 mm	200 mm	24/40	79,000
C18-10-546	L1209A27	ø50 mm	250 mm		87,300
C18-10-552	L1209A28	ø75 mm	300 mm		134,600



### High Efficiency Air Condenser 고효율 에어 냉각기

Cat. No.	Model	Length	Socket / Cone Size	EA Price ₩
C18-132-897	1218-H214	200 mm	14/40	723,400
C18-132-900	1218-H224		24/40	760,700
C18-132-903	1218-H324	350 mm	24/40	791,800
C18-132-906	1218-H424	450 mm		816,700

- 물을 사용하지 않는 냉각기
- Glass Bellow 디자인으로 넓은 냉각 표면적을 만들어줌
- Liebig Condenser와 비슷한 냉각 효율을 제공함
- Clear Visibility of the Reaction Reflux
- 5ml ~ 1L 용량의 Solvent와 같이 사용하기에 알맞음
- 진공상태에서 사용하기 적합
- Chemglass

