

Cold Finger Condenser, Cooling Rod 콜드 핑거 냉각기, LukeGL®

Cat. No.	Model	길이	NS Cone	EA Price ₩
C18-10-558	L122214	100 mm	24/40	31,100
C18-10-564	L122201	150 mm		31,100
C18-10-570	L122203	200 mm		31,100
C18-10-576	L122205	250 mm		31,100
C18-10-582	L122206	300 mm		33,700
C18-10-588	L122207	350 mm		33,700


- Borosilicate Glass 3.3
- with Cone and Water Connections
- 길이는 Joint 밑으로의 길이임
- Joint Size는 요청에 의해 변경 가능함
- LukeGL®



Micro Cold Finger Condenser 마이크로 콜드 핑거 냉각기

Cat. No.	Model	Length	Joint Length	Joint NS	EA Price ₩
C18-132-120	CG-1222-A-01	119 mm	5 mm	14/20	293,300
C18-132-122	CG-1222-A-02	127 mm	13 mm		293,300

- Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- 사용 Flask 용량 : 5 ~ 50ml
- 1 ~ 20ml 의 solvent를 Reflux 하기에 알맞음
- Hose Connection : OD 10mm
- Chemglass




Reflux Condenser, Thimble Shape 팀블형 환류 냉각기

- Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- Compact 하고 고효율의 낮은 높이를 가진 환류 냉각기
- Double Wall Thimble 모양의 냉각 라인
- Chemglass

with Glass Hose Connection, OD 10 mm ①

Cat. No.	Model	Size	Joint	EA Price ₩
C18-132-240	1217-20	w65×h225 mm	24/40	402,300
C18-132-242	1217-21	w65×h230 mm	29/42	540,500
R02-132-648	1217-25	w100×h315 mm	45/50	854,200



with Removable Hose Connection, PP Screw Cap ②

Cat. No.	Model	Size	Joint	Cap	EA Price ₩
C18-132-246	1217-HC-20	w64×h210 mm	24/40	GL14	523,100
C18-132-248	1217-HC-21	w64×h215 mm	29/42		536,500
C18-132-250	1217-HC-25	w100×h315 mm	45/50	GL18	973,600

