

## Tubing Connector

### 튜브 연결관

• Tommyheco

#### 재질 : Natural Polypropylene

- -13 ~ 90 °C, Non Autoclavable
- USP Class VI, Rohs

#### Straight Connector, 일자형 연결관, PP ①

Cat. No.	Model	Connector od (외경)	EA Price ₩
C20-140-004	PPS01	1/16" / 1.59 mm	1,500
C20-140-008	PPS02	3/32" / 2.38 mm	1,500
C20-140-012	PPS03	1/8" / 3.18 mm	1,500
C20-140-016	PPS04	5/32" / 3.97 mm	2,900
C20-140-020	PPS05	3/16" / 4.76 mm	2,900
C20-140-024	PPS06	1/4" / 6.35 mm	1,900
C20-140-028	PPS08	5/16" / 7.94 mm	1,900
C20-140-032	PPS09	3/8" / 9.53 mm	1,900
C20-140-036	PPS12	1/2" / 12.7 mm	1,900
C20-140-040	PPS16	5/8" / 15.88 mm	1,900
C20-140-044	PPS19	3/4" / 19.05 mm	3,300
C20-140-048	PPS25	1" / 25.4 mm	4,000

#### Cross Type Connector, 십자형 연결관, PP ②

Cat. No.	Model	Connector od (외경)	EA Price ₩
C20-140-052	PPC01	1/16" / 1.59 mm	3,900
C20-140-056	PPC02	3/32" / 2.38 mm	4,000
C20-140-060	PPC03	1/8" / 3.18 mm	4,000
C20-140-064	PPC06	1/4" / 6.35 mm	5,000

#### Elbow Connector, L-형 연결관, PP ③

Cat. No.	Model	Connector od (외경)	EA Price ₩
C20-140-068	PPE01	1/16" / 1.59 mm	2,800
C20-140-072	PPE02	3/32" / 2.38 mm	2,800
C20-140-076	PPE03	1/8" / 3.18 mm	3,100
C20-140-080	PPE04	5/32" / 3.97 mm	5,000
C20-140-084	PPE05	3/16" / 4.76 mm	5,000
C20-140-088	PPE06	1/4" / 6.35 mm	2,900
C20-140-092	PPE08	5/16" / 7.94 mm	3,200
C20-140-096	PPE09	3/8" / 9.53 mm	3,200
C20-140-100	PPE12	1/2" / 12.7 mm	3,200
C20-140-104	PPE16	5/8" / 15.9 mm	3,700
C20-140-108	PPE19	3/4" / 19.05 mm	3,900
C20-140-112	PPE25	1" / 25.4 mm	5,200

#### Y- Type Connector, Y-형 연결관, PP ④

Cat. No.	Model	Connector od (외경)	EA Price ₩
C20-140-116	PPY01	1/16" / 1.59 mm	3,800
C20-140-120	PPY02	3/32" / 2.38 mm	3,800
C20-140-124	PPY03	1/8" / 3.18 mm	3,800
C20-140-128	PPY04	3/16" / 4.76 mm	6,200
C20-140-132	PPY05	1/4" / 6.35 mm	4,600
C20-140-136	PPY06	3/8" / 9.53 mm	4,500
C20-140-140	PPY08	1/2" / 12.7 mm	7,500
C20-140-144	PPY09	5/8" / 15.9 mm	7,500
C20-140-146	PPY10	21.6 mm	4,700

#### 재질 : PVDF 테프론

- -62 ~ 121 °C, Autoclavable
- Ultimate Chemical Resistance
- USP Class VI, Rohs
- Certified NSF / ANSI61

#### PVDF 테프론

Cat. No.	Model	EA Price ₩
C20-140-296	PVS01	2,100
C20-140-300	PVS02	2,200
C20-140-304	PVS03	2,300
C20-140-308	PVS04	4,000
C20-140-312	PVS05	4,300
C20-140-316	PVS06	5,500
C20-140-320	PVS08	5,800
C20-140-324	PVS09	6,400
C20-140-328	PVS12	7,500
C20-140-332	PVS16	8,600
C20-140-336	PVS19	10,900
C20-140-340	PVS25	14,100

#### PVDF 테프론

Cat. No.	Model	EA Price ₩
C20-140-344	PVC01	7,000
C20-140-348	PVC02	7,000
C20-140-352	PVC03	7,000
C20-140-356	PVC06	7,000

#### PVDF 테프론

Cat. No.	Model	EA Price ₩
C20-140-360	PVE01	3,300
C20-140-364	PVE02	3,500
C20-140-368	PVE03	3,900
C20-140-372	PVE04	6,200
C20-140-376	PVE05	6,600
C20-140-380	PVE06	6,600
C20-140-384	PVE08	7,000
C20-140-388	PVE09	7,000
C20-140-392	PVE12	7,500
C20-140-396	PVE16	9,800
C20-140-400	PVE19	12,300
C20-140-404	PVE25	15,500

#### PVDF 테프론

Cat. No.	Model	EA Price ₩
C20-140-408	PVY01	4,700
C20-140-412	PVY02	4,700
C20-140-416	PVY03	4,700
C20-140-420	PVY04	7,500
C20-140-424	PVY05	8,700
C20-140-428	PVY06	22,100
C20-140-432	PVY08	25,000
C20-140-436	PVY09	32,900

